

Actim Fecal Blood

INSTRUCTIONS FOR USE

Numbers **[1]**, **[12]** refer to illustrations.

STRUCTURE OF DIPSTICK
[1] Dip area **[2]** Result area **[3]** Test line **[4]** Control line

INTENDED USE
The Actim[®] Fecal Blood Test is a visually interpreted, qualitative immunochromatographic dipstick test for detecting occult blood in feces. The Actim Fecal Blood test detects human hemoglobin in a feces suspension. The test is intended for professional use to help diagnose gastrointestinal bleeding.

KIT COMPONENTS
The components of the Actim Fecal Blood test are:

- Actim Fecal Blood dipstick, in a perforating device. The dipsticks (10) are packed in an aluminum container (30361ETAC).
- Specimen Dilution Buffer tubes (30372ETAC). This phosphate buffered solution contains bovine serum albumin (BSA), protease inhibitors and <0.1% (w/v) sodium azide (NaN3). One sampling stick is attached to the cap of each tube.

The kit 30331ETAC contains 2 aluminum tubes of dipsticks (30361ETAC) and 20 Specimen Dilution Buffer tubes (30372ETAC). The Actim Fecal Blood Buffer kit (30372ETAC) contains 20 Specimen Dilution Buffer tubes (30372ETAC).

An instructions for use leaflet is included in all product formats. Toilet paper or a clean disposable container is needed for sampling (not provided with the test).

STORAGE
Store the dipsticks and the Specimen Dilution Buffer tubes at +2...+30°C. Stored unopened, each component can be used until the expiry date marked on the component. After opening the aluminum container first time, the dipsticks can be stored for 4 months at +2...+30°C. However, the expiry date must not be exceeded.

SPECIMEN COLLECTION
The specimen is feces suspension that is prepared as soon as possible after sample collection. The original feces sample can be stored for 3 days at +2...+8°C.
Unsew the cap of the Specimen Dilution Buffer tube **[5]**. Collect the feces sample by using the sampling stick attached to the screw cap. There are two small slits in the head of each stick. Twist the stick randomly in several different places in the feces sample so that both slits become filled with feces **[6]**. In case of liquid samples, dry the stick by wiping it with a piece of toilet paper before sticking it into the sample. The slits of the stick will be filled with the liquid feces.

Put the sampling stick into the Specimen Dilution Buffer tube by pushing the stick through the cone. Close the cap tightly **[7]** and shake the tube to suspend the feces in the buffer **[8]**. The specimen suspension is now ready for testing. The suspension can be stored for one week at +2...+25°C (+2...+8°C recommended).

TEST PROCEDURE AND INTERPRETATION OF THE RESULTS
1. If stored refrigerated, allow the aluminum container of dipsticks and the buffer tube to reach room temperature before opening. Remove the required number of perforating devices (each containing a dipstick) from the container just prior to use and close the container immediately. Do not touch the yellow dip area at the lower ends of the sticks. Do not touch the part of the perforating device below the stopper, as this area will come into contact with the sample. Identifying marks may be written on the upper orange parts of the dipsticks.
2. The specimen should be at room temperature upon testing. Use the sharp tip of the perforating device to puncture the perforation area in the bottom of the buffer tube **then push the device into the tube until the stopper prevents it going further** **[9]**.
3. Slowly invert the tube (not fully up) down one and immediately place the tube on the table in an upright position **[10]**. **NOTE: It is important that the Inverting movement is not too abrupt, a peaceful inversion is vital for the correct performance of the test.**
4. The result can be interpreted as positive as soon as two blue lines become visible in the result area. The result is negative if only one blue line is visible. If the suspension has reached the lower part of the result area **[11]**, **Do not pay attention to any lines appearing later than 10 minutes.**
5. If **two blue lines, the test line and the control line**, appear, the test result is **positive**. If **one blue line, the control line**, appears, the test result is **negative**.
If the control line does not appear, the test is invalid **[12]**.

In case of a missing control line, please see the Notes section for troubleshooting.
DISPOSAL OF THE TEST DEVICE
The contact between the Specimen Dilution Buffer tube cap and the perforating device is so close that the test device can be disposed of directly after the result has been read. Do not remove the perforating device from the tube.

LIMITATIONS OF THE TEST
The tests is intended for *in vitro* diagnostic use only.
NOTES
• The test results are qualitative. No quantitative interpretation should be made based on the test results.
• Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
• Do not use a dipstick if you notice a blue coloring in the result area before testing.
• Take care that any part of the test device that will be in contact with the sample tube do not become contaminated by external blood, e.g. from a cut finger. Do not touch the dip area of the dipstick or the area of the perforating device below the stopper.
• When inverting the sample tube to allow the dipstick to absorb sample, do not invert the tube too fast, or for too long. The test will not work properly if the amount absorbed is too small or too large.
• If the control line does not appear, the inversion of the sample tube may have been performed too abruptly. Tap a new dipstick and repeat the inversion in a controlled way.
• Improper sampling may lead to false result.
• Occult blood is not uniformly distributed in the feces sample. Therefore, it is important to dip the sampling stick randomly in several places in the feces sample.
• The test line is in the lower half, the control line in the upper half of the result area of the dipstick. Appearance of a control line confirms correct performance of the test. If a control line does not appear the test is invalid, and should be repeated using another dipstick.
• If the test result cannot be interpreted clearly (e.g. if the lines are blotched or uneven) it is recommended that the test be repeated with another dipstick.
• At ten minutes the appearance of any faint-to-dark blue test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 10 minutes.

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PRINCIP OF THE TEST
Feces do not normally contain measurable amounts of blood. The presence of hemoglobin in feces implies that there may be bleeding in the gastrointestinal tract.
The test is based on immunochromatography. It involves two monoclonal antibodies to human hemoglobin. One is bound to latex-labeled particles (the detecting label). The other is immobilized on a carrier membrane to catch latex particles and indicate a positive result. When the feces suspension comes into contact with the dip area of the dipstick, the dipstick absorbs the liquid, which starts to flow up the dipstick. If the sample contains hemoglobin, it binds to the antibodies attached to the latex particles. The particles are carried by the liquid flow and, if hemoglobin is bound to them, they bind to the catching antibodies. A blue line (test line) will appear in the result area if the concentration of hemoglobin in the sample exceeds the detection limit of the test. A second blue line (control line) confirms correct performance of the test.

PERFORMANCE OF THE TEST
Analytical Sensitivity
The analytical sensitivity, detection limit, of the Actim Fecal Blood test was evaluated using samples with different concentrations of human hemoglobin on three different lots of Actim Fecal Blood test. The detection limit of the Actim Fecal Blood test is approximately 50 µg of human hemoglobin in the sample and the results remain positive at least to 500 000 µg/L. A feces suspension with 50 µg of hemoglobin/1 corresponds to 25 µg of hemoglobin/5g of feces.


Repeatability and reproducibility
The test reproducibility (repeatability) assay precision were evaluated using samples containing 0 – 2000 µg/L of human hemoglobin in feces suspensions. For repeatability the samples were tested on the same day in 10 replicates per sample on three different lots of Actim Fecal Blood test. The total amount being 180 tests in one day. For reproducibility the samples were tested on five consecutive workdays using two replicates per sample on three different lots of Actim Fecal Blood test, the total amount being 180 tests in five days. Repeatable and reproducible results were obtained.

Analytical specificity, cross-reactivity
The Actim Fecal Blood test was tested with following animal hemoglobin: Bovine, Equine, Rabbit, Goat and Dog. No cross reactivity was observed. The Actim Fecal Blood test is specific for the human hemoglobin.

Interference testing
Interference of dietary iron supplement and dietary C-vitamin supplement in the Actim Fecal Blood test was examined. The feces samples were collected before medication and during the medication. No interference of the tested substances was observed with the performance of the Actim Fecal Blood test.

Diagnostic performance
Actim Fecal Blood test was evaluated for diagnosing colorectal disease. The Ostávků klinická biochemie a 4. interní klinika VFN aL.F.UK Prague, Czech Republic. The results are shown in FIG.1.

Method comparison
The Actim Fecal Blood test was compared using the new and old Specimen Dilution Buffer tube. The agreement with the two methods was 100%.

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