EN
INSTRUMENT
MANUAL

Actim[®]
1ngeni

Copyright

This manual, as well as the hardware and software described in it, is provided under licence and may be used and/or copied only in accordance with the terms of such licence. Actim is a trademark of Actim Oy.

Notice regarding serious incidents

For a patient/user/third party in the European Union and in countries with an identical regulatory regime (IVDR, EU 2017/746/EU) if, during the use of this device or as a result of its use, a serious incident has occurred, please, report it to the manufacturer and to your national authority. The contact information for the manufacturer of this device to report a serious incident is:

Actim Oy, Kloivpellontie 3, FI-02180 Espoo, Finland Tel. +358 9 547 680
www.actimtest.com, support@actimtest.com

WARNINGS

 **Warning** Always operate the Actim 1ngeni instrument on a surface that is level, dry and not exposed to direct sunlight.

 **Warning** Do not move the instrument when conducting a test.

 **Warning** Do not touch the screen with sharp objects.

 **Warning** Do not drop the instrument.

 **Warning** Do not place objects on top of or directly against the instrument.

 **Warning** Do not operate the instrument if the power supply pack is wet.

 **Warning** Do not operate the instrument with a damaged power supply pack.

 **Warning** The Instrument Check (IC) cartridge is sensitive to light. Follow the specific storage instructions for the IC cartridge.

 **Warning** Only use approved listed peripheral accessories with the instrument.

 **Warning** Do not use this instrument in close proximity to sources of strong electromagnetic radiation (e.g. unshielded intentional RF sources), as these may interfere with its proper operation.

 **Warning** If a USB extension cable is used with the instrument, the total length should not exceed 3.0 meters.

 **Warning** The maximum current supplied by all USB peripherals should not exceed 1.1A (2 high power loads and 1 low power load).

 **Warning** This instrument is designed to operate only with the provided power supply plug pack. This module forms part of the system. Do not operate the system with a different power supply module. The correct power supply is required to maintain the safety and electromagnetic compatibility of the system.

 **Warning** Risk of electrical shock. Do not operate the instrument or the power supply plug pack if it has been opened, damaged or exposed to moisture, condensation or rain. The external power supply plug pack is sealed with no user serviceable parts. Do not operate this module with any damaged or exposed parts.

 **Warning** Do not open or attempt to repair the instrument or other accessories as there is a risk of damage to the instrument. This instrument does not contain serviceable parts and should be returned to the supplier for repair. Opening the instrument will also void the warranty.

The real time clock coin battery included in the equipment runs for the operational life and is not a user replaceable item.

 **Warning** Operate the instrument only for its intended purpose and in accordance with this user manual and warnings. The protection provided by the equipment may be impaired if the equipment is operated in a manner contradictory to the above. This instrument (including the power supply) is designed to operate within the manufacturer's specifications. Do not exceed the manufacturer's specifications during use.

 **Warning** Position the instrument with clear access to connectors. Keep the connected cables clear of work areas such that tripping or catching will not pull the instrument off of its workbench. The mains socket outlet intended to be used with the Actim 1ngeni instrument's external power pack should be located near the equipment and should be readily accessible.

 **Warning** USB and Ethernet Interfaces. If intended for connection to external equipment, please ensure that the interfaces of such equipment are separated from the mains by double reinforced insulation and present no risk of electrical shock.

 **Failure to follow these warnings will void the warranty of the instrument.**

Warnings

This device is designed for use in a professional healthcare facility environment. It is likely to perform incorrectly if used in a home healthcare environment. If it is suspected that performance is affected by electromagnetic interference, correct operation may be restored by increasing the distance between the device and the source of the interference.

Before taking Actim 1ngeni instrument into use, please read the instrument's manual and the following warnings and descriptions carefully. The Actim 1ngeni instrument is intended to provide safe and reliable operation when used in accordance to this user manual. If the instrument is used in a manner that is not specified in the user manual, the protections provided by the equipment may be impaired. The instrument is designed to operate safely under the following conditions:

ITEM	DESCRIPTION
Number of test bay slots	Configured for 1 bay.
Multiple assays	Supports multiple Actim test types and protocols.
Measurement technology	Advanced image acquisition and analysis.
Modes	Timed reading and manual reading
Illumination type	Green LED - 520 nm Red LED - 622 nm Blue LED - 470 nm
Colour Touch Screen:	3,5" diagonal capacitive TFT LCD
Communications:	Ethernet for data transfer. USB for data export, barcode reader and printer.
Data storage:	High capacity on-board storage can store up to 999 test results. History Records allow for search and retrieval. Archive or export via USB/ Ethernet.
Power:	12 V DC from external AC/DC supplied plug pack. Input: 100~240V~50/60Hz 0.4A Battery power (optional). DC Voltage fluctuation $\pm 10\%$ DC Current consumption: 0.2A DC typical @ 12V DC
Dimensions:	124 mm X 114 mm X 118 mm
Weight:	Approx. 650 g
Printer Support:	USB attached Seiko SLP-620 and SLP-650 printers.
Operating Environment:	Indoor Use (protected from water) 15°C to 30°C. 20% to 70% RH (non-condensing). 0 to 2000m altitude.
Storage Environment:	2°C to 45°C, 20% to 70% RH (non-condensing) for at least 7 days.

- Set up the instrument on a stable, level bench, in an office or laboratory environment. Position the device so that it is easy to turn power off and unplug the power supply.
- The Actim 1ngeni instrument is not intended as a hand held instrument; only operate it on a flat surface.
- Install the instrument at least 100mm from an edge.
- Install cables so as to prevent the risk of tripping or pulling that may cause instrument or personal damage.
- The instrument is a non-serviceable part, opening the instrument will void the instrument warranty.
- Ensure that ferrites are fitted to the USB peripheral accessories before operation with the instrument.
- Installation Categories (Overvoltage categories) II
- The mains socket for the AC/DC power pack should be readily accessible.

INTRODUCTION

1.1 Intended purpose

The Actim® 1ngeni is an automated rapid *in vitro* diagnostic instrument. The instrument has been designed to provide quantitative and/or qualitative results of immunochromatographic colorimetric tests measuring various analytes from patient samples to aid in diagnosis and treatment monitoring. The instrument is intended to be used by healthcare professional in clinical laboratories and near-patient testing (NPT) settings. The instrument is intended to be used with the Actim 1ngeni test kits. Target population is defined by the associated assays. Specimen type is based on the associated assays. For *in vitro* diagnostic use only.

1.2 Package components

- Actim 1ngeni instrument (Cat. No 19101AC)
- Power supply
- Adapters for international use
- Fitting ferrites for USB accessories
- Instrument Check (Cat. No. 19116IC). Instrument check is available also separately.

Materials required but not provided are listed in Actim 1ngeni test kit instructions for use.

1.3 Actim 1ngeni overview


Fig 1.
Front view of the instrument.

Fig 2.
Rear view of the instrument.

1.4 User interface

All operator interactions with the Actim 1ngeni instrument are performed using the instrument's touch screen.


1.4.1 Home screen

Displays buttons for access to the test, results and settings. The current set time and date, as well as user ID, are displayed.

1.4.2 Keyboard and numeric keypad

Onscreen keyboard and numeric keypad enable user text input into the instrument. The text can also be entered using a barcode reader.

1.4.3 Navigation tools

Returning to the home menu

Down

Log out

Back

Up

Test

Cancel

Search tool

Print

Ok

Administrator

Export

1.4.4 Menu structure

The home menu consists of three sections for the user: test (Section 5), results (Section 7) and settings (Section 3). In addition to these, the administrator (admin) will have access to the admin settings (Section 4).


Fig 3.
Menu structure

1.4.5 Power button as instrument status indicator

The different power button signals indicate the instrument's status.

Start-up: Holding down the power button for 1 second will initiate the Actim 1ngeni instrument's start up. Self test is performed in every start up assuring functionality of the instrument. Wait for the instrument self test to complete, after which the instrument's application home screen is displayed.

! Do not insert instrument check or test cartridge into the instrument during the self test.

Shutdown: Once the instrument is powered on and the instrument's application is visible, the instrument is shutdown when the power button is pressed and held for 3 seconds. The shutdown status bar is displayed and the instrument is powered off. The instrument will go into a screen saver mode when left on and not in use.

Rapid flash: Initial power connection registered.

Slow flash: Power connected. Instrument shutdown.

Solid: Instrument on.

1.5 Materials required

- Cartridge:** Used for inserting a test into the instrument. The cartridge includes a QR code with test-specific information. The cartridge is provided as part of the Actim 1ngeni test kit.
- Ethernet Cable:** Is required to connect the Actim 1ngeni to a network. Used for connectivity to Desktop Software and remote access to test results. Connect the network cable prior to powering on the instrument.
- USB Key:** Is used for importing and exporting data (not supplied with the instrument).


Fig 4.
Cartridge with QR code.

GETTING STARTED

This section assists the user and the admin user in the initial setup of the instrument.

2.1 Initial Admin user setup

At initial setup, configuration of the instrument, done by the administrator, is needed. Figure 5 shows the required steps.

Unpack the Actim 1ngeni instrument and set it up on a stable, level bench, in a clean office or lab type environment.

Configure the power supply for your region and connect the 12V power supply to the power connection of the instrument.

After the instrument is turned on, it performs automatically a self test that lasts for some minutes (see Section 6.1 for more information).

After the self test, log in as admin user:

Admin ID: "admin"

Admin password: "admin"

! Change the admin password after the first log in (see section 4.11 for more information).

Select admin icon to enter the admin settings of the instrument menu.

By using the arrow icons, scroll through the list and perform the following actions:

- Select the correct language
- Reset admin password
- Set correct time and date (Take into consideration the daylight saving times at your location).

After these basic settings define/create:

Test mode: patient test/read only (see Sections 5 and 4.3).

User list: add new users manually or by using a USB stick and define whether user ID check is on or off (see Section 4.4).

Quality test settings: define the schedule and outcome of instrument check and control test (see Section 6).

Software update and importing test types are required in case of updating the instrument (see Section 10).


Fig 5.
Admin setup
steps.

2.2 Initial user setup

After the admin user has set up the background settings, the steps shown below are defined by the user.

To start:

- Connect the 12V power supply to the power connection of the instrument.
- Log into the instrument

Enter settings by selecting the settings icon (See Section 3 for detailed instructions).


Fig 6.
Basic user
setup steps.

Define the audio settings: both touch tone and general tones used by the instrument can be adjusted from the levels 0 to 5. It is also possible to listen to the test sound.

Define the network settings: this allows you to configure the Ethernet connection of your instrument

Define the printer setup: this allows you to configure the printer connections of your instrument

Define screen brightness: adjust from 30% to 100%.

(See Section 3 for more information).

SETTINGS

The “Settings” screen is accessed from the home menu and allows a user to configure the Actim 1ngeni instrument. The different settings can be viewed using the arrow icons.

3.1 Audio

The volume can be adjusted from 0 (no sound) to 5 (highest volume) at the audio screen. The default volume is 3. It is possible to listen to the chosen sound by pressing the word “Test”.

The instrument uses sound for the following functions:

- Touch tone: A screen “click” whenever the touch screen is used
- Tones: Process tones for end of test, alerts and warnings

3.2 Network

The “Network Setup” screen allows a user to configure the Ethernet connection of the instrument.

Supported functions include:

Selecting between DHCP and Static Network

Entering IP address

Entering Subnet Mask

The default is DHCP. Switching to Static IP, IP Address and a Subnet Mask must be entered manually.

DHCP (Router): The Actim 1ngeni instrument supports the use of a typical network with DHCP enabled. Plug the instrument into a network port on the router, then power the instrument on.

Using Static: IP Addressing is a method of setting up a network where the IP address is assigned manually. The benefits of a static address is IP addresses can be assigned carefully so that each instrument has its own address—with no overlap and the IP address is fixed. This does mean however that when you connect an Actim 1ngeni instrument to a new network, you would have to select the “manual” configuration option and enter in.

3.3 Printer setup


The following printer settings can be modified:

- Select Auto Print.

3.4 Brightness

The screen brightness may be adjusted using the up and down arrows.

The default setting is 80%.



3.5 About

The Info screen displays the following information:

Serial Number, Silicon ID, MAC Address, IP Address, Application Version, Model, Baseboard Hardware Version, Core Hardware Version, Baseboard Firmware Version, Trayboard Firmware Version, Revision ID, RAM (MB), Test Types Package, Self Test Calibration Exposure, Instrument Check

! A record of these identifying fields can assist in the provision of support and assistance.

! A record of this information can be exported from the instrument using the “Export Test Results” or the “Export Summary File” and can be found on the attached USB key in the Instrument Details.csv file.

ADMIN USER SETTINGS


Admin user settings allow the admin user to configure the Actim 1ngeni instrument functions.

Enter the admin user ID.

Enter the admin user password (see Section 2.1 for initial login and Section 4.12 for instructions on changing the password). The admin password is required to access the Administrator settings. Authorised admin users will be issued the password.

The admin user sees the same Home Menu as a standard user with the addition of the admin user settings icon (see Fig 3).

Use the administrator settings screen to change or update the instrument configuration, including:

- Export
- Test Setup
- Users
- Instrument Check
- Control Test
- Self Test
- Import Test Types
- Language
- Set Time
- Set Date
- Change Admin Password

4.1 Export

The “Export” settings screen allows an admin user to export instrument data onto an attached USB key.

“**Export .csv File**” exports a .csv file containing a summary of all of the test results stored on the instrument. These results are viewable in Excel.

“**Export .TSV File**” exports a .tsv file containing a summary of all of the test results stored on the instrument. These results are viewable in Excel.

The user will be given the option to delete all of the test results after export has been performed. The exported results on the USB stick should be renamed before exporting new results.

“**Export Test Diagnostics**” enables the admin user to export data where each field value is separated from the next by a “tab” stop character (e.g. ‘2 56’). In this case, 2 and 56 will be placed in separate columns.

The user may wish to “Export Instrument Log” and provide this file to the manufacturer during service, maintenance and repair circumstances. The default setting for this function is “off”.

The Actim 1ngeni instrument holds up to 999 results. A warning will appear before the memory becomes full. The test and QC results must be exported before the memory becomes full, or memory will be overwritten.

4.2 Save Diagnostics

Turning on Save Diagnostics greatly increases the memory usage of the instrument, as it collects information in the JSON file for troubleshooting purposes. Only use as directed by instrument support personnel. The default setting is off.


4.3 Read Only

In the Read Only test mode, the test reaction occurs outside of the instrument. The user is responsible for the timing and start of the analysis. The Read Only option is visible in the test options list, if turned on by the admin user. All tests run in Read Only mode are marked in the test results. The Read Only mode is useful when tests are read in series.

4.4 Users

The Admin user has full control of the User List including following functions:

- Viewing the user list
- Adding new users
- Turning user ID check on or off (default setting is “on”)
- Exporting and importing user lists using a USB memory stick


4.5 Instrument Check

The admin user can configure the instrument check method and schedule:

The method can be configured between:

- None:** whereby the instrument check functionality is hidden from the testing and results, and no instrument check status is applied.
- Warning:** whereby a warning is displayed to the user prior to running a control test or patient test when the instrument check status is set to fail or is due.
- Lockout:** whereby the user is prevented from running a control test or patient test when the instrument check status is set to fail or is due.

The frequency at which the instrument demands the test can be chosen by the admin user (none, daily, weekly or monthly). The default frequency for the instrument check is daily, and the recommended minimum frequency is at least weekly.


(See Section 6 for more information).

4.6 Control Test

The admin user can configure the control test method and schedule:

The method can be configured between:

- None:** Where the control test functionality is hidden from the testing and results, no control test status is applied.
- Warning:** Where a warning is displayed to the user prior to running a patient test when the control status for that test type is set to fail or is due.
- Lockout:** Where the user is prevented from running a patient test when the control status for that test type is set to fail or is due.


The frequency at which the instrument demands the test can be chosen by the admin user (none, daily, weekly or monthly). The default frequency for control testing is “none”.

(See Section 6 for more information).

4.7 Self Test

The admin user can run a self test as desired and set up a schedule for automatic self test. The schedule can be chosen among the options: none, daily, weekly or monthly. After the number of days has elapsed since the last self test was run, the self test runs again.

! A scheduled self test does not interrupt testing.

(See Section 6 for more information).

4.8 Import Test Types

The instrument is supplied by the manufacturer with a default set of test types loaded onto it. If new test types are available from the manufacturer, the admin user can import test types onto the instrument. The appropriate file type needs to be loaded from a USB key.

! The USB key must be connected to the instrument before attempting the import.

If the test types file was successfully imported the instrument displays a confirming “Info” message to the user. The instrument will reboot.


Information about the test types package is viewable on the instrument’s “About” screen (see Section 3.5 for more information).

! Importing a new test types package will replace the test types package previously loaded into the device.

! The Workspace test types package must be in the USB flash memory key root directory folder for the Actim 1ngeni instrument to find the test types package file.

4.9 Language

The “Language” settings screen allows an admin user to change the language of the instrument. All text strings on the instrument will be displayed in the language selected. Text strings that are input in the test types package are independent of language settings.


The default language setting is UK English.

4.10 Set Time and Date

The Time & Date screen allows the admin user to set the correct time and date. Use the touch screen to choose between a 12 Hour or 24 Hour display.

Use the touch screen to select hour, minutes, day, month, or year, and use the up and down arrows to adjust.

Press [OK] to save the settings. The instrument will reboot in order to adopt the new date and time.

NOTE 1: The time and date is displayed to the user on the main screen. Time and date are confirmed by the user.


NOTE 2: The ACTIM 1NGENI instrument Real Time Clock (RTC) maintains the instrument set time whilst the instrument power down.

! Take into consideration the daylight saving times at your location.

4.11 Change Admin Password


The admin password can be changed by choosing “Change Admin Password” option from the “Admin Settings” menu.

Enter Current Password: The admin user must enter the current administrator password to prove they have authorisation to change the admin password.

Enter New Password: The admin user must enter a new password.

! Please ensure that the admin user password is kept in safekeeping to prevent an instrument settings lockout.

4.12 Data Connectivity

Data connectivity refers to the ability to connect to a laboratory information system (LIS) and send assay results for storage in that system.

The Actim 1ngeni Instrument communicates with an LIS Server using the HL7 standard. This standard defines TCP/IP data exchange and protocol frames used for communication between a LIS server and an instrument.

Note!. Only unidirectional communications are possible using the LIS transmission feature from the Actim 1ngeni Instrument to the LIS Server.

Packet Format

Results to be sent to the LIS are packaged as an ASCII formatted data packet. The packet is divided broadly into several categories/sections each of which contain fields appropriate for certain types of information.

As an example, a data packet may contain the following categories:

- Message Header - Information used to parse the message
- Specimen Segment - Information about the specimens tested
- Observation Request Segment - Information about the Type of test requested (1 to many)
- Observation Results Segments - Information about the result of the tests.

The categories above are defined in the HL7 standard.

Setup LIS Connection

The following section demonstrates the steps for Setting up the Actim 1ngeni for LIS communications.



Step 1. Select LIS Setup from the Data Connectivity menu

The first step is to Set the LIS Server IP address and Port Number settings in the Actim 1ngeni to establish communications on the same network.

Step 2. Set LIS Server IP Address

Setting up communications to the LIS Server requires knowledge of the server credentials. Contact your IT provider to acquire the IP Address & Port Number details of the LIS Server to be used.

- Select "IP Address"
- Enter the IP Address of the LIS server into the field provided, with each unique ID number separated by a full stop '.' E.g. "XXX.XXX.XXX.XXX"

Step 3. Set LIS Server Port Number

Setting up communications to the LIS Server requires knowledge of the server credentials. Contact your IT provider to acquire the IP Address & Port Number details of the LIS Server to be used.

- Select the “Port Number” option.
- Enter the Port Number of the LIS Server into the field provided.

The default value is provided is 51112. However, this number will be specific to the LIS Server you are using. The number must lie between the following range: 49152 - 65535

Step 4. Test LIS Server Connection

To Test the connection between the Actim 1ngeni & the LIS Server, run the ‘Test Connection’ option:

- Select “Test Connection”
- One of three states will be displayed:
 - Pending: Test Connection not yet run
 - Success: Actim 1ngeni successfully connects to LIS
 - Failed: Actim 1ngeni failed connection with LIS

Touch the OK icon, to apply confirm the applied IP & Port Number settings of the LIS Server.

Use the ‘Clear Settings’ option, to quickly clear all configured options within the LIS Setup menu.

Check Connectivity Status

- Connectivity Status

The user may access the Status screen to view transmission data statistics with the LIS Server.

To do this, navigate to the Data Connectivity menu under the Admin settings, then select the ‘Connectivity Status’ option

- LIS Status

Under the LIS Status menu, the user can review the transmission stats with the LIS Server.

The categories shown are:

- Unsent Results (results not yet sent to LIS)
- Total Results (total number of results in instrument memory)

LIS Status	
Unsent Results:	500
Total Results:	654
Test Connection:	Pending
Send All Unsent	

- Test connection

The first user option shown is a repeat command also found under the LIS Setup menu.

The user can run the ‘Test Connection’ option from the LIS Status screen as well for convenience.

- Send All Unsent

The 'Send All Unsent' option performs an immediate transmission of all unsent test results to the LIS.

Refer to Section 1.4.4 Toolbar Indicators, to understand the Transmission Status shown by the Actim 1ngeni Instrument.

When sending the results, the unsent result count will count down to 0 as each previously unsent result is sent to the LIS server.


During transmission of the results, a stop button (X) will appear on the bottom right corner of the screen. This will allow the use to stop the current transmission and leave the remaining results as unsent.

Auto Sending

The user can set the Actim 1ngeni Instrument to automatically send results to the LIS Server at regular periods. When the period elapses, all of the unsent results will be sent to the LIS Server.

Auto send can be configured by selecting on the "Auto Send" option from the Data Connectivity menu within the Admin Settings.

Home > Admin Settings > Data Connectivity > Auto Send.

Toggle the Auto Send ON / OFF button , to enable the Auto Send Period
Adjust the Auto Send period by using the Up/Down arrows.

The available periods are:


- 5 Minutes
- 15 Minutes
- 30 Minutes
- 1 Hour
- 1 Day

Press the OK icon to confirm the applied settings.

Clear Sent Status

The user can clear the sent status of the number of transmissions sent to the LIS Server.

This may be useful if the user exports all test results off the instrument and clears the memory and wishes to track one for one any new test results.

Send Single Result to LIS Server

There are two methods to send a single result to the LIS Server:

- At the end of running a live test,
- Via the Results menu


At the end of running a test select the Export icon  to open the result export options.

When viewing the Result page select the Export icon  to open the result export options.

Select the "Send to LIS" button to transmit the currently viewed result to the LIS server; regardless of its previous sent/unsent status.

NOTE: LIS must be configured before the user is able to Export Result to LIS

NOTE: The export button is disabled from result screens if LIS is not configured and no USB is detected

RUNNING A TEST

It is essential that the on-screen prompts are followed throughout the testing process.

5.1 Patient test

Select a new test by pressing the test icon on the home menu. Enter a Patient ID (1-20 characters) using the onscreen keyboard or a barcode scanner (see Section 9 for more information).

The Actim 1ngeni instrument displays the test number and the patient ID and asks the user to insert the test cartridge. For sampling and activating the test, refer to the Actim 1ngeni test instructions for use document. The test begins when cartridge insertion is detected. Selecting back cancels the test and no data is saved. The test type and user ID are saved.




The test screen will display:

- Lot number for the test
- Expiry date for the test
- Test number
- Patient ID
- Test type
- User ID
- Test time and date

Remaining time of incubation

! An expired test cannot run on the instrument.

The test result will either be positive or negative (see Section 11 For error messages). The cartridge can be detached from the instrument, when the result window is visible.

The patient result window includes:

- Test number
- Patient ID
- Test type
- LIS status
- Time and date of the Test
- Lot number of the Test
- Expiry date
- User ID
- Control line result
- Result

The result can be printed or exported to an attached USB Key (see Section 7 for more information) or LIS (see Section 4 for more information)

5.2 Patient Test (Read Only)

In the read only test mode, the test reaction occurs outside of the instrument. The user is responsible for the timing and start of the analysis.

Select a new test by pressing the test icon. Enter a Patient ID (1-20 characters) using the onscreen keyboard or a barcode scanner (see Section 9 for more information).

The Actim 1ngeni instrument displays the test number and the patient ID and asks the user to insert the test cartridge. For sampling and activating the test, refer to the Actim 1ngeni test Instructions for Use document. The test result is read when cartridge insertion is detected. Selecting back cancels the test and no data is saved.


! An expired test cannot run on the instrument.

The test result shall either be a positive or negative result (see Section 11 for error messages).

The patient result window includes:

- Test Number
- Patient ID
- Time and Date of Test
- Test Type (Read Only) mode
- Lot Number for the test
- Expiry Date
- User ID
- Result


The result can be printed or exported to an attached USB Key (see Section 7 for more information) or LIS (see Section 4 for more information)).

QUALITY CONTROL

6.1 Self test

The Self Test is an automated instrument test that occurs automatically during the instrument start up, or as scheduled by the admin user. It can also be run by the admin user at any time. Running self test independently from the start up and determining the schedule for an automated test can be done using the admin settings.

! Do not insert instrument check or test cartridge into the instrument during the self test.

After the self test, the following status screens will apply:

Pass: In the start up self test case, the instrument enters the login screen.

Warning: "Self Test Warning!" is displayed. In the start up self test case: the user confirms and advances to login screen and testing is not locked out. When requested or scheduled by the admin user: the user confirms and returns to previous application screen. Testing is not locked out.

Fault: "Self Test Fault! Testing has been locked out" is displayed. In the start up self test case: the user confirms and advances login screen. When requested or scheduled by the admin user: the user confirms and returns to previous application screen. Patient and control testing is locked out.

The current self test result is available in the instrument information screen and the last passed self test is shown in the printed report output.

Each failed Self Test is assigned a level of severity:

Fault: the Self Test Status is set to fail and the testing is locked out.


Warning: the Self Test Status is set to fail but the testing is not locked out.

Critical Failure: the Self Test Status is set to fail and the instrument reboots.

6.2 Instrument check

The Instrument Check (IC) function provides the capability to periodically check the functions of the instrument. This check has the advantage that it is independent of any assay and assay controls and specifically checks the instrument reading capability using an external cartridge and printed strip standard.

Determining the outcome and frequency of the instrument check can be done at the admin settings.

Instrument check can be run by pressing the test icon and choosing instrument check.

➔ INSTRUMENT CHECK


The instrument displays the test number and asks the user to insert the test. The test begins when cartridge insertion is detected.

! Use a special cartridge containing the IC test strip, designed only for the instrument check.

The test result is either pass or fail. Failed results do not display a test outcome, instead display “Error: {0}” as the test result. The error code can be looked up in the Error section of this user manual. In case of a failed result, it is recommended to turn the instrument power off and on again, and rerun the instrument check using an unused IC cartridge.

The Instrument Check result window includes:

- Test Number
- Time and Date of Test
- Lot Number for the IC cartridge
- Expiry Date for the IC cartridge
- User ID
- Result
- LIS status
- Control line result

The result can be printed or exported to a USB key (see Section 7 for more information).

6.3 Liquid controls

➔ CONTROL TEST

Actim liquid control kits are available separately. When the instrument is used for the first time, running successful positive and negative QC runs is recommended. If the control testing is set to lockout by the admin user, the user is required to run both a successful positive and negative QC test before a patient test is allowed. Liquid controls may also be tested at any time in order to conform to local requirements.

The liquid control test can also be performed in read only mode. The control status for individual test types can be viewed through the control test menu.

Enter a unique control ID using the onscreen keyboard or a barcode reader to input the text field. Press the OK icon to save the text input (1-20 characters).

The instrument displays the test number and control ID and asks the user to insert the test. The test result is read when cartridge insertion is detected. Selecting “back” cancels the test and does not save any test data.

The user must select the associated control test with a known outcome that they wish to run for the inserted test. Selecting “Cancel” cancels the test and does not save any test data.

! The test will timeout if a selection is not made before the timer ends.

The instrument displays the lot number and expiry date for the test inserted into the instrument in addition to the test number and control ID.

The test time is taken from the cassette insert event.

The test type, user ID and remaining time of incubation are shown on the display.


! An expired test cannot run on the instrument.

The test result are either pass or fail. Failed results do not display a test outcome, instead display “Error: {0}” as the test result. The error code can be looked up in the Error section of this user manual.

The control result window includes:

- Test Number
- Control ID
- Time and Date of Test
- Control Type
- Lot Number for test used
- Expiry Date for the test used
- User ID
- Result
- LIS status

The result can be printed or exported to an attached USB Key (see Section 7 for more information) or LIS (see Section 4 for more information).

6.4 Control status

The Control Status for each test type (for example, PROM or Partus) is indicated by an icon.

! Failed quality controlling results may have an effect on the test results.

Fail

Fail is shown when one or more of the associated control tests have failed.

Due

Due is shown when one or more of the associated control tests are due.

Pass

Pass is shown when all of the associated control tests have passed.

RESULTS

7.1 Patient Test Results

The instrument automatically saves test results to a location in the instrument's memory or to LIS if configured (See Section 4 for more information).



The results menu is entered with the result icon:

From "Patient Test", a result list of all patient tests saved in the instrument memory is shown. The following information is displayed:

- Patient ID
- Test Date and Time

The results summary screen can be navigated through using the up and down navigation arrows.


7.2 Search for results

The user can search for particular patient test results from the result list by entering one of the following search terms:

- Test Number
- Patient ID
- Time and Date of Test
- Test Type
- Lot Number
- Expiry Date
- User ID

7.3 Printing and exporting results

  **PATIENT TEST** (e.g.)  (export) /  (print)

The user can select the result from the results summary screen to see the individual result. The patient result window includes:

- Test number
- Patient ID
- Test type
- LIS status
- Test time and date
- Lot number
- Expiry date
- User ID
- Control line result
- Final Test result

The result can be printed or exported to an attached USB Key (see Section 7 for more information) or LIS (see Section 4 for more information).

! Remember to plug in the desired accessory or USB key before pressing the export or print icon.

7.4 Control Test Results

  **CONTROL TEST RESULTS**

If Control Tests are selected from the results menu then a summary list of all of the Control Tests saved in the instrument memory will be displayed with the following information:

- Control Test Type
- Test Date and Time

The results summary screen can be navigated through using the up and down navigation arrows.


The user can filter the summary list by entering a search term (see Section 7.2 for more information).

The user can select the result from the results summary screen to see the individual result.

The Control result displays:

- Test Number
- Control Lot ID
- Time and Date of test
- Test Type
- LIS status
- Lot Number (if available)
- Expiry Date (if available)
- User ID
- Control line result
- Result

(See Section 7.3 for more information on printing and exporting the results).

7.5 Instrument Check Results


> INSTRUMENT CHECK

If Instrument Check is selected from the results menu then a summary list of all of the Instrument Checks saved in the instrument memory will be displayed with the following information:

- Instrument Check result
- Test Date and Time

The results summary screen can be navigated through using the up and down navigation arrows.


The user can filter the summary list by entering a search term (see Section 7.2 for more information).

The user can select the result from the results summary screen to see the individual result.

The result displays:

- Test Number
- Time and Date of test
- Test Type (Read Only) mode
- Lot Number
- Expiry Date
- User ID
- Result

(See Section 7.3 for more information on printing and exporting the results).

CLEANING, MAINTENANCE AND SERVICE

⚠ WARNING: The isopropyl alcohol used in this procedure is flammable.

Ensure that the instrument is not powered.
Do not use isopropyl alcohol within 3m of an open flame or source of ignition.
Avoid skin contact.

⚠ WARNING: Instrument may be contaminated.

Avoid skin contact.
Wash hands with hand wash after completing decontamination

All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

Suggested materials:

Gloves:	Disposable laboratory gloves
Wipes:	Lint free wipes
Swabs:	Foam tipped swab
Isopropyl Alcohol:	99% isopropyl alcohol in a spray dispenser
Hand wash:	Disinfectant hand wash

The Actim 1ngeni instrument can be cleaned using a lint free wipe, damp with isopropyl alcohol (IPA). Do not use free liquids to clean the instrument.

Dip the lint free wipe into the isopropyl alcohol and allow any excess fluid to flow off. If any lint or dust remains on the wipe, dispose of the wipe and use a new one.

Inspect:	Inspect for damage or visible contamination.
Dispose:	Dispose of any materials left on the instrument such as test parts.
Wipe surfaces:	Wipe all surfaces of the instrument with wipes moist with isopropyl alcohol.
Dispose:	Dispose of all used materials and gloves.
Wash hands:	Wash hands with the disinfectant hand wash.

Service

Actim 1ngeni is free from regular maintenance. In case of instrument malfunction or repair requirement, please contact your local supplier. Delete all patient results from the instrument before shipping it to service.

Disposal

Actim 1ngeni instrument is a low-voltage electronic device. A used Actim 1ngeni instrument must be handled as potentially biohazardous waste. The instrument should be disposed of as Used Medical Equipment according to national and local law. Delete all patient data before disposing of the instrument. The packaging materials are recyclable.

Revision history

Revision history can be found from www.actimtest.com

OPTIONAL EQUIPMENT

Optional equipment that can be used with the Actim 1ngeni instrument can be purchased separately from their supplier:

- USB key
- Label printer, the Seiko SLP620 and SLP650
- Barcode scanner, the Datalogic QuickScan Barcode Wand QD2430

In addition to these, it is possible to use an external power pack (Tekkeon MP3450).

Three ferrite cores are supplied with the Actim 1ngeni instrument. The ferrite must be fit into the following accessories prior to the use of the instrument:

Place the USB peripheral cable inside the ferrite core; the distance from the ferrite core to the USB connector base should be $29\pm 2\text{mm}$

Lock the cable into place by pushing the cable down.

Close the ferrite core once the cable is in the correct position.

SOFTWARE UPDATE

 **Caution:** Stored test data may be at risk of being deleted during the software update process! It is highly recommended that the test result stored on the Actim 1ngeni instrument is archived to an external formatted, USB Key prior to performing a software update

 **Caution:** DO NOT REMOVE THE USB KEY DURING THE SOFTWARE UPDATE. The USB Flash memory instrument can be removed from the instrument once the software update is complete.

Software update:

1. Export all of the results data from the Actim 1ngeni instrument.
2. Place the update file in the main directory of a clean USB key.
3. Insert the USB key containing the update into the instrument USB port.
4. Power off the instrument.
5. Power on the instrument and follow the on-screen update instruction. Touch to start.
6. When update is complete, remove the USB key. Touch to continue is displayed.
7. The instrument will reboot after the update is completed.
8. Observe the version number displayed during start up to confirm that the software update was applied.

Check also the about screen on the instrument to confirm update version number.

 **The first time a self test is run following a software update, a fault or warning result may be given. Please make special note of the information given on the self test screen as to what other functions may need to be conducted following the software update. Please run the self test for a second time from the admin menu.**

TROUBLESHOOTING

This section provides troubleshooting steps for:

1. In-Test Error Dialogs
2. Warning Dialogs
3. Error Dialogs
4. Information Dialogs

Once the steps listed below are executed, if the error or warning persists, please contact your local distributor.

Warning: The software uses a standard warning screen. A warning screen will be displayed where a user has made a selection that is not reversible and a confirmation is required.

Error: If the instrument receives a request or performs an action outside of the normal operating parameters for that user/test/function, an error message will be displayed, explaining what error has occurred and requiring the user to confirm that the error has been noted before the instrument returns to normal operation. In some cases it may cancel a current test and/or require an instrument reboot.

Information: An information screen supplies the user with important information. There are no "OK" or "Cancel" inputs required from the user, so once the user has finished reading the message "OK" will exit the screen.

In-Test Error Dialogs

CODE	DESCRIPTION	ACTION
		<p>Possible cause: the strip is not present in the strip carrier. Check that the strip is in the cartridge/strip carrier. Check for contaminants. Attempt to rerun test</p>
Error: 2	<p>Could Not Identify Cartridge This can occur only if the overall image average grey level is below a defined limit.</p>	<p>Possible cause: the camera or LED's have failed and the image is black Run a self test on the instrument. If the self test fails, contact the supplier</p> <p>Possible cause: exposure calibration is not correct – re-execute the exposure calibration. Run a self test on the instrument. If the instrument needs to be calibrated then contact the supplier.</p>
Error: 5	<p>Could Not Locate Strip The instrument could not locate the cartridge/strip carrier in the acquired image.</p>	<p>Check that strip is correctly inserted in the cartridge/strip carrier. Change it to check for contaminants on the strip carrier/ cartridge fiducial marks. Attempt to rerun test Try another strip carrier/cartridge. If the error continues, contact the supplier. "Acquire image" of the instrument internals and check the image quality is adequate.</p>
Error 20	<p>Could not normalise</p>	<p>The signal from the strip could not be 'normalised' during the test. Typically, this occurs when the cartridge/carrier is located in the wrong position in the image. Ensure that your cartridge is being inserted correctly and repeat the test. If this error occurs repeatedly contact the test package developer or the manufacturer.</p>
Error: 36	<p>Cartridge/Carrier Could Not Be Correctly Found - Scale The instrument fiducial features on the cartridge/strip carrier found but the scale is out of range.</p>	<p>Check that the correct cartridge/strip carrier is being used with the instrument. Check that the cartridge/strip carrier was inserted correctly into the instrument. Check for contaminants. Check that the selected test matches the cartridge/carrier used. Attempt to rerun test using a new cartridge/strip carrier. If the error continues, contact the supplier.</p>
Error: 37	<p>Cartridge/Carrier Could Not Be Correctly Found - Position The instrument fiducial features on the cartridge/strip carrier found but position is out of range.</p>	<p>Check that the correct cartridge/strip carrier is being used with the instrument. Check that cartridge/strip carrier was inserted correctly into the instrument. Check for contaminants. Check that the selected test matches the cartridge/carrier used. Attempt to rerun test using a new cartridge/strip carrier. If the error continues, contact the supplier.</p>
Error: 38	<p>Failed to locate the strip cartridge/ carrier The instrument fiducial features on the cartridge/strip could not be found. The test was unable to proceed.</p>	<p>Check that the correct cartridge/strip carrier is being used with the instrument. Check that cartridge/strip carrier was inserted correctly into the instrument. Check for contaminants. Check that the selected test matches the cartridge/carrier used. Attempt to rerun test using a new cartridge/strip carrier. If the error continues, contact the supplier.</p>

In-Test Error Dialogs

Error: 41	<p>Cartridge/Carrier Could Not Be Correctly Found - Rotation The analysis software was able to locate the cartridge/strip carrier; however the rotation angle found was outside of the acceptable limits.</p>	<p>Check that the correct cartridge/strip carrier is being used with the instrument. Check that the cartridge/strip carrier was inserted correctly into the instrument. Check for contaminants. Check that the selected test matches the cartridge/carrier used. Attempt to rerun test using a new cartridge/strip carrier. If the error continues, contact the supplier.</p>
Error: 43	<p>Could not locate the control line, Multiple Candidates The instrument was unable to determine the location of the control line with sufficient confidence. The algorithm found more than one control line. This error may occur when no control line is present on the strip.</p>	<p>Check that strip is correctly inserted in the cartridge/strip carrier. Check for control line being present. Check the position of the lines on the strip to ensure that it is within manufacturing specifications/tolerances. Check that a test line is not within the search region for the control line. Check for contaminants on the strip. Repeat testing with new dipstick. If the error continues, contact the supplier. The Supplier will review the test types package.</p>
Error: 44	<p>Strip Analyser Could Not Locate Control Line, No Line Found The instrument has found no control lines within the control line search region.</p>	<p>Check that the strip is correctly inserted in the cartridge/strip carrier. Check for the control line. Check the position of the lines on the strip to ensure that it is within manufacturing specifications/tolerances. Check for contaminants on the strip. Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package.</p>
Error: 45	<p>Strip Analyser Could Not Locate Line, Line Width Validation Failed A line has been found, but the line width score is outside of the acceptable limit defined in the test type. Either the test strip/cartridge is invalid, or the test type is incorrectly defined.</p>	<p>Check that strip is correctly inserted in the cartridge/strip carrier. Check the test line. Check for contaminants on the strip. Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package.</p>
Error: 46	<p>Strip Analyser Could Not Locate Line, Line Peak Validation Failed A line has been found, but the peak score is outside of the acceptable limit defined in the test type. Either the test strip/cartridge is invalid, or the test type is incorrectly defined.</p>	<p>Check that strip is correctly inserted in the cartridge/strip carrier. Check the test line. Check for contaminants on the strip. Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package.</p>
Error: 47	<p>Strip Analyser Could Not Locate Line, Line Area Validation Failed A line has been found, but the area score is outside of the acceptable limit defined in the test type. Either the test strip/cartridge is invalid, or the test type is incorrectly defined.</p>	<p>Check that strip is correctly inserted in the cartridge/strip carrier. Check the test line. Check for contaminants on the strip. Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package.</p>
Error: 48	<p>Strip Analyser Could Not Locate Line, Line Position Validation Failed A line has been found, but the position score is outside of the acceptable limit defined in the test type. Either the test strip/cartridge is invalid, or the test type is incorrectly defined.</p>	<p>Check that strip is correctly inserted in the cartridge/strip carrier. Check the test line. Check for contaminants on the strip. Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package.</p>

In-Test Error Dialogs

Error: 49	<p>Focus Out Of Range The focus of the camera is outside of limits. This error will only occur if the test has the 'Run Focus Check' option enabled in kinetic.</p>	<p>Re-power the instrument. Attempt to run a self test. Supplier to acquire image and ensure that the internal validation marks are clearly visible and not obscured by debris or other materials.</p>
Error: 50	Analyser Exception	<p>A general analysis error. Re-power the instrument. Attempt to run a self test. If error continues contact supplier. "Acquire image" of the instrument internals and check the image quality is adequate</p>
Error: 52	No Image Acquired	<p>Camera failed to return an image. Re-power the instrument. Attempt to run a self test. If error continues contact supplier. "Acquire image" of the instrument internals and check the image quality is adequate.</p>
Error: 64	<p>Could Not Make Decision, Unknown Exception The decision module in the test type encounters an error.</p>	<p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. Other Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package.</p>
Error: 65	<p>Could Not Make Decision, Audio Tone Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The Audio Tone is not being set by the algorithm and is required.</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the audio tone is set correctly by the algorithm.</p>
Error: 66	<p>Could Not Make Decision, Detailed Message Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Detailed Message' Fields are not set by the algorithm.</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the Detailed Message fields are set correctly by the algorithm.</p>
Error: 67	<p>Could Not Make Decision, Icon Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Icon' fields are not set by the algorithm. (RESULT_DECISION_ICON_X)</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the result Icons are set correctly by the algorithm.</p>

In-Test Error Dialogs

Error: 68	<p>Could Not Make Decision, Message Not Set</p> <p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Message' fields are not set by the algorithm. (RESULT_DECISION_MESSAGE_X)</p>	<p>Possible causes are:</p> <p>The test package is not correctly installed – try importing again.</p> <p>An incorrect test package is being used with the software version installed on the instrument. Check test package release details.</p> <p>Contact the supplier.</p> <p>The Supplier will review the test types package to ensure that the result messages are set correctly by the algorithm.</p>
Error: 69	<p>Could Not Make Decision, Title Not Set</p> <p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Title' fields are not set by the algorithm (RESULT_DECISION_TITLE_X)</p>	<p>Possible causes are:</p> <p>The test package is not correctly installed – try importing again.</p> <p>An incorrect test package is being used with the software version installed on the instrument. Check test package release details.</p> <p>Contact the supplier.</p> <p>The Supplier will review the test types package to ensure that the result 'Title' fields are set correctly by the algorithm</p>
Error: 70	<p>Could Not Make Decision, Set UI Type Not In Allowed List</p> <p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Type' field is not set to a valid value by the algorithm.</p>	<p>Possible causes are:</p> <p>The test package is not correctly installed – try importing again.</p> <p>An incorrect test package is being used with the software version installed on the instrument. Check test package release details.</p> <p>Contact the supplier.</p> <p>The Supplier will review the test types package to ensure that the UI type for the result is set correctly by the algorithm.</p>
Error: 71	<p>Could Not Make Decision, UI Type Not Set</p> <p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Type' field is not set by the algorithm.</p>	<p>Possible causes are:</p> <p>The test package is not correctly installed – try importing again.</p> <p>An incorrect test package is being used with the software version installed on the instrument. Check test package release details.</p> <p>Contact the supplier.</p> <p>The Supplier will review the test types package to ensure that the UI Type is set correctly by the algorithm.</p>
Error: 72	<p>Could Not Make Decision, Unknown Type</p> <p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Type' field is set to an invalid type for the application.</p>	<p>Possible causes are:</p> <p>The test package is not correctly installed – try importing again.</p> <p>An incorrect test package is being used with the software version installed on the instrument. Check test package release details.</p> <p>Contact the supplier.</p> <p>The Supplier will review the test types package to ensure that the UI Type is set correctly by the algorithm.</p>
Error: 73	<p>Could Not Make Decision, Valid Flag Not Set</p> <p>This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Valid' field is not set to true by the algorithm.</p>	<p>Possible causes are:</p> <p>The test package is not correctly installed – try importing again.</p> <p>An incorrect test package is being used with the software version installed on the instrument. Check test package release details.</p> <p>Contact the supplier.</p> <p>The Supplier will review the test types package to ensure that the UI Type is set correctly by the algorithm.</p>

In-Test Error Dialogs

<p>Error: 74</p>	<p>Could Not Make Decision, Ratio Title Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Ratio Title' field is not set by the algorithm.</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the required number of ratio output fields are defined correctly.</p>
<p>Error: 75</p>	<p>Could Not Make Decision, Ratio Output Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Ratio Output' fields are not set by the algorithm</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the required ratio outputs are set correctly.</p>
<p>Error: 76</p>	<p>Could Not Make Decision, QC Result Not Set Correctly This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'QC Result' output is not set by the algorithm correctly.</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the QC result flag is set correctly.</p>
<p>Error: 80</p>	<p>Could Not Make Decision, Quant Title Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Quantitative Title' output fields are not set correctly.</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the quantitative output fields are set correctly.</p>
<p>Error: 81</p>	<p>Could Not Make Decision, Quant Output Not Set This is most likely due to an error in the logic of the decision algorithm definition in the test package. The required 'Quantitative Output' fields are not set correctly.</p>	<p>Possible causes are: The test package is not correctly installed – try importing again. An incorrect test package is being used with the software version installed on the instrument. Check test package release details. Contact the supplier. The Supplier will review the test types package to ensure that the quantitative output fields are set correctly.</p>
<p>Error: 82</p>	<p>Strip Analyser Could Not Locate Control Line, No Peak Candidates The strip analyser could not locate any trace of a control line in the expected control line region.</p>	<p>Check the control line. Check for contaminants on the strip. Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package to ensure that the control line search region is sufficiently wide enough to enclose the control line.</p>
<p>Error: 83</p>	<p>Strip Analyser Could Not Locate Test Line, No Peak Candidates The strip analyser could not locate any trace of a test line in the expected test line region.</p>	<p>Attempt to rerun test If the error continues, contact the supplier. The Supplier will review the test types package, ensure that the test line search region is sufficiently wide enough to enclose the test lines.</p>

In-Test Error Dialogs

Error: 86	Normalisation Exception	The normalisation failed to complete. Re-attempt normalisation. Contact the Manufacturer.
Error: 96	Could Not Complete Exposure Calibration, Not Scale Calibration	Could not complete exposure calibration, the exposure calibration algorithm is not able to complete the primary scale calibration. Check that the strip is in the cartridge/strip carrier. Check that the cartridge has fiducials markings that are visible and not damaged. Check for contaminants. Attempt to re-run test.
Error: 145	Strip Analyser Scale Rotate Factor Out Of Range	Could not analyse, the image scale and rotation is out of range. The instrument requires re-calibration
Error: 200	Could Not Locate Cartridge, No Edges Could not locate the cartridge window in the expected location.	Check that cartridge is inserted correctly into the slot. Check for damage or contaminants on cartridge. Attempt to re-run test
Error: 201	Could Not Locate Cartridge, No Label	Could not locate the cartridge printed label in the correct locations. Check that cartridge is inserted correctly into the slot. Check for damage or contaminants on cartridge. Attempt to re-run test.
Error: 300	Image Level Out Of Range	The camera returned a dark or light image that is out of range. Re-power the instrument. Attempt to run a self test. If error continues contact supplier. "Acquire image" of the instrument internals and check the image quality is adequate
Error: 301	Voltage Out Of Range	Per-test self test failed. The system voltages are out of range. Re-power the instrument. Attempt to run a self test. If error continues contact supplier. Ensure correct power supply is being used
Error: 302	Temperature Out Of Range	Per-test self test failed. The instrument temperature is out of range. Re-power the instrument. Attempt to run a self test. The ambient air temperature is out of the acceptable range. If error continues contact supplier.
Error: 310	Consumable Removed	The consumable was removed during image acquisition.
Error: 320	Barcode Not Found	The barcode could not be found in the image. Ensure that the correct cartridge is being used. Check that cartridge is inserted correctly into the slot. Check for damage or contaminants on cartridge Attempt to re-run test.
Error: 321	Barcode Has Expired	The barcode on the cartridge has expired. Discard cartridge
Error: 322	Barcode Does Not Match Selected Test	The barcode on the cartridge does not match the selected test type. Attempt to re-run test with the correct test type selected.

In-Test Error Dialogs

Error: 323	Barcode Check Digit Not Correct	The barcode check digit is incorrectly formatted. Check that cartridge is inserted correctly into the slot. Check for damage or contaminants on cartridge. Attempt to re-run test. Discard cartridge if error continues.
Error: 324	Barcode In Incorrect Location	The barcode was read, but the cartridge is not in the correct location for the test. Check that cartridge is inserted correctly into the slot. Attempt to re-run test.
Error: 325	Barcode Format Error	The format of the barcode is incorrect. Check that cartridge is inserted correctly into the slot. Change for damage or contaminants on cartridge. Attempt to re-run test. Discard cartridge if error continues.
Error: 326	Barcode Found But Damaged	The barcode has been found but is damaged. Check that cartridge is inserted correctly into the slot. Change for damage or contaminants on cartridge. Attempt to re-run test. Discard cartridge if error continues.
Error: 145	Strip Analyser Scale Rotate Factor Out Of Range	Could not analyse, the image scale and rotation is out of range. The instrument requires re-calibration.

Warning Dialogs

0001	Are you sure you want to delete user {0}? This cannot be undone.	The instrument requires confirmation that the user does intend to delete the user ID. This step once confirmed cannot be undone.
0002	Delete lot data for Test ID: {0} Lot #: {0} Expiry: {0}	The instrument requires confirmation that the user does intend to delete the selected lot data. This step once confirmed cannot be undone.
0003	Are you sure you wish to change the Admin password?	The instrument requires confirmation that the admin user wishes to continue.
0004	Do you wish to change the date to: {0}? This may affect the status of the instrument schedulers.	The instrument requires confirmation that the user wishes to change the instrument date. This is required because changing the instrument date will affect the self test, instrument check and control test schedules.
0005	Do you wish to change the time to: {0}? This may affect the status of the instrument schedulers.	The instrument requires confirmation that the user wishes to change the instrument time. This is required because changing the instrument time will affect the self test, instrument check and control test schedules.
0006	Import test types package: {0} Version = {1}? Import replaces existing test types.	A instrument warning informing the user that the Test Types file they are attempting to install on the Actim 1ngeni instrument will replace the current test type set.
0007	The control status for this test has failed or is overdue! Do you wish to proceed?	The instrument is warning the user that the instrument either; failed its last Control test for that test type or it is due to for another Control test. The user may proceed with the test but the test result may be compromised.
0008	Instrument check has failed or is overdue! Do you wish to proceed?	The instrument is warning the user that the instrument either; failed its last instrument check test or it is due to for another instrument check test. The user may proceed with the test but the test result may be compromised.

Warning Dialogs

0009	Reset instrument settings to factory default? All current settings and data will be lost.	The instrument requires confirmation from the user that they wish to remove all user entered data and patient results from the instrument by returning the instrument to factory default.
0010	Instrument memory is almost full! Please export and delete result data.	This dialog notifies the user that the Actim 1ngeni instrument memory is almost full and needs to clear in the near future. Attach a USB flash memory key to the instrument and perform Test Result export and confirm the request to wipe instrument memory. This dialog will appear until the memory is cleared or the memory is full.
0011	Instrument memory is full! {0} is deleted if you wish to proceed.	This dialog notifies the user that the Actim 1ngeni instrument memory is full and needs to clear before further testing can be performed. If you proceed a result will be deleted from the test memory to make room for the new test result. This cannot be undone. Attach a USB flash memory key to the instrument and perform Test Result export and confirm the request to wipe instrument memory.
0013	Delete all test results from the instrument memory? This cannot be undone.	The instrument requests the user to confirm that the user would like to delete all test results from the instrument. This cannot be undone
0014	Cancel changes? Changes made will be lost.	This is a instrument warning, advising the users that any changes made on the current GUI screen will be lost by navigating away without saving changes.
0015	Cancel test? Test data will be lost.	This is a instrument warning, advising the users that cancelling the current test will lose the current test data. This cannot be undone.
0017	Import list of {0} users? Imported user list will replace existing user list. This cannot be undone.	The instrument informs the user that importing a new user list will replace the user list on the instrument. Please make sure this is what you want to do before proceeding with the import request.
0022	A remote connection is still active! Are you sure you wish to cancel?	If the Actim 1ngeni instrument has been used in Connected mode with the Actim 1ngeni Kinetic Software then it may find the connection still open. To close the connection, select "OK"; to keep it open select "Cancel".
0023	Save diagnostics substantially increases memory use! Do you wish to proceed?	The Save diagnostics function on the instrument uses more memory than a standard test result. Only use this function for troubleshooting purposes. The instrument warns you before you turn this function on.
0025	Leaving this page will reset its values. Are you sure you want to proceed?	This is an Actim Oy factory user item only.
0026	Are you sure you wish to proceed?	This is a secondary confirmation requesting the user confirm that they do indeed wish to return the instrument to factory default.
0029	Cancelling will lose the normalisation progress made! Are you sure you want to cancel?	This is an instrument warning, advising the users that cancelling the current normalisation progress will lose the current data. This cannot be undone.
0032	Reset device settings to user default? All current settings and data will be lost	The instrument requires confirmation that the admin user wishes to reset the instrument to a user default state. This cannot be undone.

Warning Dialogs

0035	Are you sure you wish to change the password?	The instrument requires confirmation that user wishes to change the user password.
0036	Are you sure you want to clear the LIS settings?	The instrument requires confirmation that the user wishes to clear the current LIS settings
0038	Import Test Types Package: {0} Version = {1} Import shall replace existing Test Types.	The instrument informs the user that importing a new test type package will replace the current test type package on the instrument. Please make sure this is what you want to do before proceeding with the import request.

Errors Dialogs

0512	A critical error has occurred! Please refer to user manual. Reboot required.	The Actim 1ngeni instrument has suffered a critical error. The instrument will not boot up. Arrange for the return of the instrument.
0513	RTC failure! Please refer to user manual.	The Actim 1ngeni instrument real time clock battery has failed. Arrange for the return of the instrument.
0514	User ID not recognised! Please try again.	The user ID entered does not match an ID entered in the instrument user ID list. Please re-attempt user ID entry. If you have forgotten your user ID, please contact the administrator.
0515	User ID input invalid! Entry should be between 1 and 20 characters.	The text input by the user does not meet the requirements of being between 1 and 20 alphanumeric characters.
0516	Invalid input! Entry should be between 1 and 20 characters.	The text input by the user does not meet the requirements of being between 1 and 20 alphanumeric characters.
0517	User ID already exists! Please enter a different ID.	The user ID entered already exists on the instrument, please either: Enter a different User ID Delete current user ID Edit current user ID ! The following user ID are not available on the instrument "admin", "factory" or "Actim Oy".
0518	Lot Expired.	If an expired Lot is detected then the instrument will not permit the activation of the Test Lot.
0519	Starting remote mode failed! Please try again.	If the user attempts to establish a remote connection with the Actim 1ngeni instrument connected via network to the kinetic software and the attempt fails, it is suggested that repeated attempts are made once the user establishes that the network and physical hardware setup is configured correctly.
0520	Barcode does not match test selection.	If an internally read barcode test selection does not match a known test type on the instrument the instrument will not permit the test to proceed.
0522	Passwords do not match! Please re-enter password.	The password entered does not match the password saved on the instrument. Please try to attempt entering the password again. If you have forgotten your password, please contact the supplier.
0523	Import of user list failed! Please try again.	The import of the user list failed. Please check that the USB key is correctly connected to the instrument and re-attempt import.

Errors Dialogs

0524	Export of user list failed! Please try again.	The export of the user list failed. Please check that the USB key is correctly connected to the instrument and re-attempt export.
0525	Barcode could not be read	The instrument was unable to read an internal barcode from a carrier/cartridge. The instrument will not permit the test to proceed.
0526	Import failed! Multiple test type packages detected. Please check USB key and try again.	The instrument has detected that there is more than one test types package available for import on the attached USB key. Remove one of the test types packages from the USB key then re-attempt import.
0527	No test type package found! Please check USB key contents and try again.	If no Test types are loaded onto the instrument then a Test Type Package must be imported onto the instrument from an attached USB Flash memory key.
0528	A maximum of 50 test types can be imported onto the instrument! Please edit the test types package	The test types package the user is attempting to import is too large. Please contact the supplier. The supplier must ensure that test types package provided to the end user contains less than 50 tests.
0529	Test types package import failed! Please try again.	The instrument was unable to import the Test Types from an attached USB Flash memory key. Please ensure that the file is placed correctly in the main file directory and the file name correct. Ensure that only a single test type package is located on the USB key.
0530	SD card not found! The instrument will reboot. Please refer to user manual.	The Actim 1ngeni instrument cannot find the external SD memory card. The instrument will not boot up. Arrange for the return of the instrument.
0532	Control status has failed! Please update the control status.	The control status for the selected test type has failed. Please run a new control test for that test type to update the control status to pass.
0533	Exposure calibration failed! Please try again.	The attempt made to calibrate the Actim 1ngeni instrument has failed. Please re-attempt calibration. If this does not correct the issue then contact the manufacturer.
0534	Normalisation failed! Please try again.	The attempt made to normalise the Actim 1ngeni instrument has failed. Please re-attempt normalisation. If the problem persists, replace the normalisation strip. If this does not correct the issue then contact the manufacturer.
0535	Instrument check has failed! Please run a new instrument check.	The user may not* be able to run a test until the instrument check test has passed. Run a new instrument check to update the instrument check status to Pass. *Instrument check functionality is configurable and the settings may vary.
0536	Self test has failed! Testing is locked out. Please refer to user manual.	The instrument self test has failed and testing has been locked out. Run a new self test to confirm, then refer to the recorded failure as to what has failed.
0539	{0} instrument exposure calibration has not been completed. Please refer to user manual.	{0} displays either Visual or Fluorescent. The instrument has not been exposure calibrated yet. Please contact the administrator to complete the exposure calibration process.
0546	Input invalid. Input must be between 1 and 20 characters.	The text input by the user does not meet the requirements of being between 1 and 20 alphanumeric characters.

Errors Dialogs

0547	Print set to {1} failed! Please check printer status.	Ensure the correct printer has been connected (SLP620). Reinsert printer and re-power. If this does not correct the issue then contact Actim Oy.
0548	No USB instrument found! Please check USB instrument connection and try again.	If the user attempts to perform a task on the instrument that required a USB Flash memory key attached, the instrument will look for the attached USB instrument. If the instrument cannot be found then an error message is displayed. Check that the USB Flash memory key is correctly attached and reattempt task.
0553	Image acquisition failed! Please try again.	If you attempt to run a test the instrument cannot take an image then the image acquisition will have failed. Please attempt to run a self test to confirm that the instrument is operating correctly.
0554	Network settings not applied! Please try again.	The system was unable to apply the selected network settings. Check network connection and try again. If the problem persists, restart the unit and try again.
0557	Test expired! Please discard the test.	If a cartridge/strip carrier is inserted into the instrument with an expired test, then the instrument will not permit the test to proceed.
0559	Control status is overdue! Please update the control status.	The control status for the selected test type is overdue. Please run a new control test for that test type to update the control status to pass.
0560	Instrument check is overdue! Please run a new instrument check.	The instrument check test status is overdue. To reset the instrument check status to pass, please run a new instrument check test.
0561	Test type not available! Please discard the test.	If the user attempts to run a test that is not available on the instrument, the instrument shall notify that the test type is not available. Please discard the test.
0562	Export failed! Please try again.	The instrument was unable to export to an attached USB Flash memory key. This could due to the following reasons: a) A USB Flash memory key was not correctly inserted into the instrument's USB serial port at the point of export. b) The USB Flash memory key was not formatted correctly and cannot be recognised by the Actim 1ngeni instrument. See the USB requirements.
0563	Results data has been corrupted! Please refer to user manual.	A results corruption has occurred, please attempt to export the results. Please contact your supplier.
0564	The test has timed out! Please discard the test.	The next test step was not completed in the allotted time. The test has timed out. Please discard the test.
0565	Password not recognised! Please try again.	The password entered is not recognised by the instrument, please enter the correct password. If you have forgotten the admin password, please contact your supplier.
0566	IP address {0} is invalid! Please enter a valid IP address.	An incorrect or invalid IP address has been entered. Please check what was entered then try again.
0567	Subnet mask {0} is invalid! Please enter a valid subnet mask.	An incorrect or invalid subnet mask has been entered. Please check what was entered then try again.

Errors Dialogs

0569	Printer not found! Please check printer connection and try again.	If the user attempts to print a test report before setting up the instrument printer connections, this error will be displayed. a) Try connecting the USB report printer to the instrument. If the printer is not found after connecting, try to reboot the instrument.
0570	Instrument is not normalised! Testing is locked out. Please refer to user manual.	If the instrument is not successfully normalised prior to testing, the test results may be incorrect. This instrument requires normalisation to be complete before the Actim 1ngeni instrument can perform testing.
0572	Maximum of 99 users reached! Please delete an existing user ID before adding a new user.	The instrument is notifying the user that the user list is full. No more users can be added until some are cleared from the instrument memory. Please delete a user to be able to add a new user.
0591	Cannot proceed!	User is unable to proceed with the test selection category as there are no test types available. If this occurs contact supplier as the test category may not be part of the supplied test type package file.
0593	Invalid Password. Passwords must be between 8 and 64 characters in length	An incorrect or invalid password has been entered. Please check what was entered then try again.
0594	Return to user default request failed! Please try again.	The instrument has not successfully been returned to user default. Re-attempt the request. If request continues to fail, run a self test. Contact supplier.
0595	Port Number Invalid.	Please enter a valid Port number. An incorrect or invalid port number has been entered. Please check what was entered then try again.
0597	Could not obtain IP Address. Please check network connection and try again.	The IP Address was unable to be obtained from the network. Ensure the instrument is connected to a valid network and try again.
0600	An error has occurred! Changes were not saved.	An error occurred writing the changes to the instrument SD card. Reboot instrument and try again. If error continues to occur, run a self test. Contact supplier.
0765	Invalid barcode found.	The barcode scanned is not recognized by the instrument. Please try again or use the keyboard to enter the details.

Information Dialog

0256	User deletion completed successfully!	The delete user function has completed successfully.
0257	Lot data deleted successfully.	The delete Lot data function has completed successfully.
0258	{0} users imported successfully!	User import function has completed successfully.
0259	Exported {0} users successfully to USB instrument!	User export function has completed successfully.

Information Dialog


0260	No user ID list found! Please check USB key and try again.	If the user is attempting to import a user list and the instrument cannot find the user list on the attached USB key then the user should check the USB key contacts to ensure that the file is in the correct location and in the correct format.
0261	No user ID list found on instrument! Please enter users.	If there are no user ID entered on the instrument, the instrument shall inform the user. The admin user should add users to the instrument then reattempt request.
0262	Test types imported successfully! Instrument shall reboot.	The instrument has successfully imported Test Types. The instrument informs the user that an instrument reboot is required prior to testing with the newly imported Test Types.
0264	Calibration completed successfully!	Calibration attempt on Actim 1ngeni instrument has been completed successfully.
0265	Normalisation completed successfully!	Normalisation attempt on Actim 1ngeni instrument has been completed successfully.
0266	Factory default settings restored successfully! Device shall reboot.	On the factory default settings have been successfully applied then the instrument must reboot.
0268	Report has been successfully sent to printer: {0}	A dialog to inform the user that the instrument successfully sent the report to the printer. If the report does not print, the issue is most likely to be with the printer.
0269	The current operation was cancelled successfully! Press OK to continue.	If the user cancels an operation, this dialog informs the user that the operation was successfully cancelled.
0270	Exported successfully to USB instrument!	The instrument successfully exported to an attached USB flash memory drive.
0272	No test results found in instrument memory!	An attempt was made to export test results when the memory was empty. Run a test and retry the export function.
0273	Image: {0} acquired successfully!	The instrument has successfully acquired an image and saved it to the attached USB key.
0276	Result deletion completed successfully!	Results deletion function has completed successfully. All results have been deleted.
0277	Changes made were saved.	Changes made were saved successfully.
0278	User ID successfully changed!	User ID was successfully changed.
0279	Password changed successfully!	Admin password was successfully changed.
0280	No results found containing {0}! Please enter different search term.	If the text in the search field entered by a user in the Test Results search function return no results then an information dialog box informs the user.
0281	Instrument time has not been set! Please set the instrument time.	Attempt to reset time.
0288	User default settings restored successfully! Device shall reboot	Upon successful restoration of the user default settings, the instrument must reboot.
0289	New user created successfully.	New user was successfully created.

Information Dialog

0291	Changes made were saved. Device will now reboot	Changes made were saved successfully. The device will reboot (when this screen is acknowledged) to adopt the changes.
0337	Time and date set successfully. Instrument will reboot.	Time and date were successfully created and instrument must reboot.
0510	Lot successfully created.	New test lot was successfully activated

SYMBOLS AND LABELS

LABEL	DESCRIPTION
	Caution, consult documents
	Consult instructions for use
	CE mark
	<i>In vitro</i> diagnostic medical instrument
	Device for near-patient testing
	Manufacturer
	Serial number
	Catalogue number
	Temperature limitation
	Keep dry
	Fragile, handle carefully
	DC
	Waste electrical and electronic equipment

 **Actim Oy**
 Klovinpellontie 3, FI-02180 Espoo, Finland
 Tel. +358 9 547 680
 www.actimtest.com
 support@actimtest.com

Actim is a trademark of Actim Oy.

WARRANTY

The Actim 1ngeni instrument is guaranteed against defects in materials and workmanship for a period of one (1) year. For specific warranty information, please contact the Actim Oy. The following defects are specifically excluded from the warranty:

Defects caused by improper operation or by the improper packaging of returned goods.

Repair or modifications done by anyone other than Actim Oy.

Materials not specified by Actim Oy.

Deliberate or accidental misuse.

Damage caused by disaster.


Damage due to use of improper solvent or sample.

The warranty does not apply to fuses.

For enquiry or request for repair service, contact Actim Oy after confirming the serial number of your instrument.

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Actim är ett registrerat varumärke av Actim Oy.
Reparation eller ändringar som utförs av någon annan än Actim Oy.
Användning av material som inte anges av Actim Oy.
Vid frågor eller begäran om reparation, kontrollera instrumentets serienummer och kontakta Actim Oy.