Actim® Partus Controls

EN

EXTERNAL CONTROLS FOR ACTIM® PARTUS TEST, INSTRUCTIONS FOR USE

INTENDED USE

Actim® Partus Controls are intended to be used with the Actim® Partus test (31921ETAC) for external quality control. The controls may also be used to demonstrate negative results and weak and strong positive results.

KIT COMPONENTS

Actim Partus Controls kit (31900ETAC) consists of one of each of the following components (please see Warnings and precautions section):

- Actim Partus Negative Control (31980ETAC)
- Actim Partus Positive Control, Low (31981ETAC)
- Actim Partus Positive Control, High (31982ETAC)
- Actim Reconstitution Solution (30889ETAC)

Instructions for use

A pipette capable to measure 500 µl of solution is needed for accurate dispensing of Reconstitution Solution (not provided with the kit).

STORAGE

Store the components at +2...+25 °C. Stored unopened, each component can be used until the expiry date marked on the component. After reconstitution, the controls should be used within 24 hours if they are kept at room temperature (below +30 °C) or within 7 days if they are kept at +2...+8 °C. Alternatively, the reconstituted controls can also be stored frozen at -20 °C for 24 months. Do not freeze after thawing.

CONTROL MATERIAL

The Positive Controls consist of purified human phosphorylated insulin-like growth factor binding protein-1 (phIGFBP-1) in a buffered protein solution with BSA and Thimerosal as preservative. The Negative Control consists of the same matrix without added antigen.

LIMITATIONS OF THE TEST

- The controls are to be used with the Actim Partus test (manufactured by Actim Oy) only.
- The controls are intended for professional in vitro use only

TEST PROCEDURE AND INTERPRETATION OF THE RESULTS

Reconstitution:

1. Unscrew the cap and open the vials by carefully removing the grey stopper.
2. Reconstitute the controls by adding 500 µl of Reconstitution Solution to the vials. When pipetting, avoid the risk of carryover by changing the pipette or its tip after each addition.
3. Close the stopper and allow the vial to reconstitute for at least 15 minutes.
4. Mix the solutions.

Testing:

1. The reconstituted controls are tested by dipping the dipsticks directly into the control vials. Each dipstick requires approximately 150 µl of sample. Before dipping, ensure that the volume of the control in the vial is sufficient.
2. Testing is performed by carefully following the instructions for use of the Actim Partus test.

The test results are interpreted as instructed in the paragraph "Test procedure and interpretation of results" of the Actim Partus test:

- The negative control gives a negative test result.
- The low positive control gives a weak positive result.
- The high positive control gives a strong positive result.

WARNINGS AND PRECAUTIONS

- Before use, ensure that the package and vials are intact. Do not use a control if the vial is damaged.
- The lyophilized controls and the reconstitution solution contain <0,01% reaction mass of: 5-chloro-2-methyl-2H-isotiazol-3-one and 2-methyl-2H-isotiazol-3-one (3:1) [EC No: 611-341-5] (Skin Sens 1, Aquatic Chronic 3). Reconstituted controls contain <0,02% reaction mass of: 5-chloro-2-methyl-2H-isotiazol-3-one and 2-methyl-2H-isotiazol-3-one (3:1) [EC No: 611-341-5] (Skin Sens 1, Aquatic Chronic 3). May cause an allergic skin reaction (H317). Harmful to aquatic life (H402). Harmful to aquatic life with long lasting effects (H412). Avoid breathing mist/ vapour (P261). Contaminated work clothing should not be allowed out of the workplace (P272). Wear eye protection/protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents/container in accordance with local regulation (P501).
- Do not use the controls before they have been completely reconstituted.
- When testing, the controls should be at room temperature to ensure optimal test results with Actim Partus test.
- Carefully follow the instructions for use of the Actim Partus test and Actim Partus Controls to ensure correct results.
- If the results of the controls are not appropriate, repeat the procedure with another dipstick.
- The frequency of quality control testing is to be determined by each laboratory independently, and will depend on factors such as patient sample frequency and/or local regulations.
- The phIGFBP-1 antigen in the Positive Controls has been shown to be negative for HBVDNA, HIV-RNA and HCV-RNA. Nevertheless, such tests are unable to prove the complete absence of viruses, and the controls should therefore be handled taking appropriate precautions.
- Like all biological materials, the controls should be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

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