
## INSTRUCTIONS FOR USE


Actim<sup>®</sup>PROM  
1ngeni

EN	Instructions for use	3	NL	Gebruiksaanwijzing	53
CS	Návod k použití	9	NO	Bruksanvisning	59
DA	Brugervejledning	15	PL	Instrukcja użytkowania	65
DE	Gebrauchsanweisung	21	PT	Instruções de utilização	72
ES	Instrucciones de uso	27	RO	Instrucțiuni de utilizare	78
FI	Käyttöohje	34	RU	Инструкция по применению	85
FR	Mode d'emploi	40	SV	Bruksanvisning	92
IT	Istruzioni per l'uso	46			

## STRUCTURE OF THE DIPSTICK

## TEST PROCEDURE AND RESULTS

## STRUCTURE OF THE CARTRIDGE

# INSTRUCTIONS FOR USE

The numbers and letters refer to the pictures on the inside cover

## STRUCTURE OF THE DIPSTICK

[A] Dip area [B] Result area [C] Test line [D] Control line [E] Handling area

## STRUCTURE OF THE CARTRIDGE

[F] QR code

## INTENDED USE

Actim<sup>®</sup> PROM 1ngeni is a qualitative immunochromatographic rapid test for detection of amniotic fluid in cervicovaginal secretions during pregnancy. The test is analyzed with the automated Actim<sup>®</sup> 1ngeni instrument. Actim PROM 1ngeni detects IGFBP-1 (Insulin-like growth factor-binding protein 1), which is a major protein in amniotic fluid and a marker of the amniotic fluid in cervicovaginal secretions. The test is intended for professional use in clinical laboratories and near-patient testing (NPT) settings to help diagnose the rupture of fetal membranes (ROM) during pregnancy when a rupture of fetal membranes is suspected. **For *in vitro* diagnostic use only.**

## BACKGROUND

Premature rupture of fetal membranes (PROM) is defined as a rupture of fetal membranes before onset of the labor. PROM is associated with increased risk of intra-uterine infection and thus it increases the risk of both maternal and perinatal morbidity and mortality. If the rupture of membranes occurs before 37 weeks of gestation, it is referred to as preterm PROM (PPROM)<sup>1</sup>. PROM occurs in approximately 5-10% of all pregnancies and PPROM in approximately 1-3 %. Management of a patient with PROM or PPROM requires evaluation of the relative risks of preterm birth and prolongation of gestation. The Actim PROM test is intended to assist the healthcare professional in making an informed decision about how to manage the patient.

The concentration of IGFBP-1 in amniotic fluid is 100–1,000 times higher than in maternal serum<sup>2</sup>. Seminal fluid and urine do not contain a measurable amount of IGFBP-1, and it is not present in the cervicovaginal secretions if the fetal

membranes are intact<sup>2</sup>. After the rupture of fetal membranes, amniotic fluid with a high concentration of IGFBP-1 mixes with cervicovaginal fluids<sup>2</sup>. In the Actim PROM 1ngeni test, a specimen is taken with a sterile polyester swab from cervicovaginal fluid, and the specimen is extracted into the extraction buffer. The presence of IGFBP-1 in the extracted specimen is detected using a dipstick.

## PRINCIPLE OF THE TEST

The test is based on immunochromatography. It involves two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles (the detecting antibody). The other is immobilized on a carrier membrane (the capture antibody) to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick is placed in an extracted specimen, the dipstick absorbs liquid, which starts to flow up the dipstick. If the specimen contains IGFBP-1, it binds to the antibody bound to latex particles. The particles are carried by the liquid flow and, if IGFBP-1 is bound to them, they bind to the capture antibody. Actim 1ngeni instrument will recognize a blue line (test line) if the concentration of IGFBP-1 in the sample exceeds the cut-off of the test. A second blue line detected by the instrument, the control line, confirms successful completion of the test.

## KIT COMPONENTS

Actim PROM 1ngeni kit contains 10 test packs (Cat. No 30831RETAC) with the instructions for use.

The components for performing one test are packed in one Actim PROM 1ngeni test pack (Cat. No 30821RETAC):

- **[SWAB]** 1 sterile polyester swab for specimen collection (EO).
- **[EXTR] [BUF]** 1 tube of extraction buffer (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives (please see Warnings and precautions section).
- **[DIP]** 1 dipstick in a sealed aluminum foil pouch with desiccant.
- **[CART]** 1 cartridge for reading the result with the Actim 1ngeni instrument. The QR code on the cartridge includes the dipstick lot number, expiry date, and lot specific calibration data.

Materials needed but not provided: Timer

## STORAGE

Store the test kit at +2...+25 °C. When stored unopened, each component can be used until the expiry date marked on both the test pack and the kit box. The test kit and the test packs can also be stored for 2 months at +25...+30 °C, as long as the expiration date is not exceeded. All reagents are ready to use. Let the reagents reach room temperature before use, if stored at +2...+8 °C. Use the dipsticks shortly after their removal from the aluminum foil pouch. Do not freeze.

Before performing a test, read and follow the test instructions for use.

## SPECIMEN COLLECTION AND EXTRACTION

If the kit components are stored refrigerated, allow the aluminum foil pouch and the extraction buffer to reach room temperature (+15...+25 °C) before use. The specimen is a cervicovaginal secretion that is extracted into the provided extraction buffer. A specimen is collected using a sterile polyester swab (provided in the kit). The specimen should be collected prior to performing digital examination and/or transvaginal ultrasound. Take care not to touch anything with the swab before taking the specimen.

1. Separate the labia and carefully insert the tip of the swab into the vagina toward the posterior fornix until resistance is met. Alternatively, the specimen can be taken from the posterior fornix during a sterile speculum examination. The swab should be held **in the posterior fornix for 10–15 seconds** to allow it to absorb the cervicovaginal secretion.
2. Open the extraction buffer tube and extract the specimen immediately from the swab by swirling the swab vigorously **in the extraction buffer for 10–15 seconds**. Press the swab against the wall of the extraction buffer tube to remove any remaining liquid from the swab. Discard the swab after extraction. **NOTE!** Do not leave the swab in the tube.

## SPECIMEN STORAGE

Specimens should be tested as soon as possible after they have been extracted, and no later than 4 hours following the specimen extraction. They should be stored at +2...+25 °C. If a specimen cannot be tested within this time, it should be frozen. The frozen specimens are stable for 24 months (at or below -20 °C). After thawing, the specimens should be mixed and tested as described below.

## TEST PROCEDURE – PATIENT TEST

In this procedure, test reaction occurs inside the instrument and results are read automatically at 5 minutes. Insertion of the cartridge into the Actim 1ngeni automatically initiates the Actim PROM 1ngeni test procedure.

3. Switch the Actim 1ngeni instrument on (if power off). After the automatic self-test, log onto the instrument by entering your user ID either manually or with a barcode reader. Select Test icon.
4. Select Patient Test icon. Enter patient ID either manually or with a barcode reader.
5. When the Actim 1ngeni screen shows “Activate and insert test”, open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch. Place the yellow dip area into the extracted specimen.
6. Hold the dipstick in the tube **until you see the liquid front enter the result area**.
7. Immediately remove the dipstick from the buffer, place it in the cartridge (reading area forward and colored side facing up) and insert the cartridge inside the Actim 1ngeni (QR code forward and facing up) until cartridge clicks into place. Follow and respond to any information appearing the screen.

Use only the cartridge included in the test pack. Ensure that the dipstick is firmly inside the cartridge. Ensure that the timer appears on the screen.

Proceed to section “Results”.

## TEST PROCEDURE – PATIENT TEST, READ ONLY

In the read only-procedure, the test reaction occurs outside Actim 1ngeni. The user is responsible for the timing and for starting the analysis. Insertion of the cartridge into Actim 1ngeni automatically initiates the reading of the test result. As the test must be read at 5 minutes, ensure that Actim 1ngeni is available for analysis before activating the test.

- Switch the Actim 1ngeni instrument on (if power off). After the automatic self-test, log onto the instrument by entering your user ID either manually or with a barcode reader. Select Test icon.

- Select Patient Test, Read Only icon. Enter patient ID either manually or with a barcode reader.
- When the Actim 1ngeni screen shows “Activate and insert test”, open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch. Place the yellow dip area into the extracted specimen.
- Hold the dipstick in the tube until you see the liquid front enter the result area.
- Immediately remove the dipstick from the buffer, place it in the cartridge (reading area forward and colored side facing up) and place the cartridge in horizontal position and start your external timer for 5 minutes. After 5 minutes, insert the cartridge immediately inside the Actim 1ngeni (QR code forward and facing up) until cartridge clicks into place. Follow and respond to any information appearing the screen, until the “Analyzing” screen appears.

Use only the cartridge included in the test pack. Ensure that the dipstick is firmly inside the cartridge.

Proceed to section “Results”.

## RESULTS

The Actim 1ngeni instrument automatically reports and saves the result. Do not visually interpret the result under any circumstances before or after analyzing the test because the test is optimized for the Actim 1ngeni instrument.

8. When the test result is visible on the screen, remove and discard the dipstick and the cartridge. The result can be either printed or exported into a USB stick or stored in Actim 1ngeni for later use. Refer to the instrument manual for details. Log off when leaving the instrument.

Positive and negative results with additional test-related information are shown on the screen. If problems have occurred during the test process, an error is reported, and no result is shown.

## ANALYZING THE CONTROL

Please refer the instructions for use provided with the control package. Analyze the control samples according to the instructions and choose “Control Test” on the display of the Actim 1ngeni instrument. The result will be stored as a control measurement.

## QUALITY CONTROL

Regular use of Actim® PROM Controls (Cat. No. 30800ETAC) is recommended. Actim PROM Controls kit includes 3 freeze-dried controls (a negative, a low positive, and a high positive control) and one reconstitution solution. The controls are provided separately.

Please consult Actim PROM Controls instructions for use and Actim 1ngeni manual for further information and instructions on how to use and analyze the controls.

Quality control shall be performed according to the applicable national or local regulations or accreditation requirements and laboratory standard Quality Controls Procedures. It is recommended to perform a control test at least with new reagent lot and shipment.

## LIMITATIONS OF THE TEST

- The test is intended for *in vitro* diagnostic use only.
- For professional use only.
- No quantitative interpretation should be made based on the test results.
- If rupture of fetal membranes has occurred **but the leakage of amniotic fluid has ceased** more than 12 hours before the specimen is taken, IGFBP-1 may have been degraded by proteases in the vagina and the test may give a negative result.
- A positive Actim PROM test result, although detecting the presence of amniotic fluid in the specimen, does not locate the site of the rupture.
- A negative test result is an indication of the present condition and cannot be used to predict the forthcoming.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.

## WARNINGS AND PRECAUTIONS

- Do not eat, drink or smoke in the area where the kit and specimens are handled.
- Wear protective personal equipment when handling and performing a test. Wash hands after specimen and test handling.
- The extraction buffer contains <0.01% reaction mass of: 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) [EC No: 611-341-5] (Skin Sens 1, Aquatic Chronic 3). May cause an allergic skin reaction(H317). Harmful to aquatic life (H402). Harmful to aquatic life with long lasting effects (H412). Avoid breathing mist/vapour (P261). Contaminated work clothing should not be allowed out of the workplace (P272). Wear eye protection/protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents/container in accordance with local regulation (P501).
- Do not use the product after the expiry date marked on the test pack or kit box. Do not exceed the stated stability periods for opened components.
- Do not mix the components from different tests and lot numbers. Components are disposable, do not reuse them.
- The test requires about 150 µl of extracted specimen to ensure proper performance of the test.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick, if you notice a blue coloring in the result area before testing.
- Do not use the dipstick, if its aluminum foil pouch or the seals of the pouch are not intact. Do not use the swab if the package is damaged.
- Use only the swab that is provided with the kit.
- Improper sampling and/or testing procedure may lead to false results.

- Care must be taken when placing the dipstick in the extraction buffer tube. The upper part of the dipstick must stay dry.
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly, it is recommended that the test be repeated.
- Do not pay attention to any lines appearing after 5 minutes.
- All biological specimens and materials must be treated as potentially hazardous and disposed of in accordance with local authority guidelines.
- Temperatures above or below the recommendations can cause inaccurate results.
- The availability of the workflow and read only modes is dependent on the administrator settings.
- A test must always be read in read only mode in case the test reaction has occurred outside the instrument.
- Do not insert anything other than the cartridges provided with Actim 1ngeni tests into the instrument.
- Proper timing of interpretation the dipstick is essential. If the interpretation of the dipstick is delayed, the test should be repeated.
- Ensure that the dipstick is firmly inside the cartridge and does not slip out when placing it into the instrument.
- If the procedural control line is not detected by the Actim 1ngeni instrument, the result will be declared indeterminate, and the test should be repeated.
- The result must not be visually inspected after or before analyzing it with the instrument.
- Do not re-analyse a test with the instrument.
- For additional information, please refer to the Actim 1ngeni instrument manual.
- In case of serious incident, please report it to the manufacturer or its representative and/or national authority.

## PERFORMANCE OF THE TEST

### Analytical Sensitivity

The analytical sensitivity of Actim PROM 1ngeni was evaluated using samples with different concentrations of IGFBP-1. The samples were tested using three replicates. The cut-off of Actim PROM 1ngeni test is approximately 25 µg/l in extracted specimen.

### Measuring range

The measuring range of the Actim PROM is approximately 25-200000 µg/l in extracted specimen.

### Analytical Specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP proteins at concentrations ranging from 10-5000 µg/l of each protein in extracted specimen using three different lots of the Actim PROM test. No cross-reactivity was seen using human IGFBP-2, -3, -4, -5 and -6 proteins. The Actim PROM test is specific to human IGFBP-1.

### Reproducibility

Reproducibility was assessed using samples containing 0 – 50 µg/l of IGFBP-1. Two operators performed the testing on different days. Three lots of the test were used; each lot read on a different instrument using the two test modes. The specimens were tested using ten replicates for a total of 120 tests per concentration. Negative samples were all identified negative and positive samples were all identified positive. One hundred percent (100 %) agreement with expected test results was obtained.

### Method comparison

Performance of Actim PROM 1ngeni was compared to Actim PROM using spiked cervicovaginal samples containing different concentration levels of IGFBP-1. The testing was performed by five operators using five different instruments. Overall agreement between the test formats was 100 %. These results confirm that the performance of Actim PROM 1ngeni test is similar to Actim PROM test (Table. 1)

TABLE 1.

	Actim PROM positive	Actim PROM negative	Total
Actim PROM 1ngeni positive	50	0	50
Actim PROM 1ngeni negative	0	50	50
Total	50	50	100

### Diagnostic performance

The Actim PROM test has been evaluated in several clinical studies. A selection of typical studies is shown in Table 2.

TABLE 2.

	Sensitivity (%)	Specificity (%)	Positive predictive value (PPV) (%)	Negative predictive value (NPV) (%)
Rutanen 1996 <sup>3</sup>	100	95	93	100
Ragosch 1996 <sup>4</sup>	100	83	83	100
Jain 1998 <sup>5</sup>	100	89	76	100
Kubota 1998 <sup>6</sup>	95	93	90	96
Ergemoglu 2004 <sup>7</sup>	97	97	97	97

### Interference testing

Take care not to touch anything with the swab before taking the specimen. Do not contaminate the swab or cervicovaginal specimen with lubricants or creams, as they may physically interfere with absorption of the specimen onto the swab and/or affect the test performance. Large amounts of meconium causing green staining of the amniotic fluid may interfere with the detection of IGFBP-1 with Actim PROM.

The following substances, conditions and micro-organisms were tested with Actim PROM test and were found not to affect Actim PROM test performance, when tested at the concentrations shown in Table 3.

TABLE 3.

Substance/Condition/ Micro-organism	Tested concentration
Whole blood	100% v/v
Semen	100% v/v
Urine	100% v/v
pH 3.5-8.5	N/A
Surgical lubricant water-based (Aquasoniq 100, Optilube, Surgilube, HR Lubricating Jelly)	50% w/v
Skin oil (Ceridal Lipolotion)	50% v/v
Pevaryl (active ingredient: econatzol.nitras)	30 mg/ml
Gyno-Trosyd (tioconazol)	20 mg/ml
Flagyl (metronidazole)	100 mg/ml
Canesten (clotrimazol) capsule	40 mg/ml
Personal lubricant water-based (K-Y Jelly, Klick)	25% w/v
Personal lubricant silicone-based (Magic Glide)	3% w/v
Baby oil (Natusan)	50% v/v
Baby powder (Natusan)	50% w/v
Feminine deodorant (Vagi-Gard)	50% w/v
Vaginal film (VCF Dissolving Vaginal Cleansing Film)	0.1% w/v
Vaginal gel (Refresh)	25% w/v
Shower and bath products (Rexona, Dove, Palmolive, Adidas, Daily)	0.1% w/v
Candida albicans	$11.2 * 10^8$ CFU/ml
Gardnerella vaginalis	$8.6 * 10^8$ CFU/ml
Neisseria gonorrhoea	$10.6 * 10^8$ CFU/ml
Chlamydia trachomatis	$5.0 * 10^3$ TCID <sub>50</sub>
HSV-1	$5.0 * 10^3$ TCID <sub>50</sub>
HSV-2	$5.0 * 10^3$ TCID <sub>50</sub>

## TRACEABILITY

The Actim PROM 1ngeni test is traceable to Behrings Institut's PP12 (IGFBP-1) preparation 307/323. This IGFBP-1 was purified from soluble extracts of human term placenta or decidua<sup>8</sup> to prepare reference standards for a quantitative in-house method, IGFBP-1 IEMA to which the Actim PROM 1ngeni result is traceable to.

## DISPOSAL

Dispose of the contents according the national and local law.

All patient specimens and used swabs, cartridges, dipsticks and extraction buffer tubes should be handled and disposed of as potentially infectious material.

Materials used in Actim PROM kit and its components:

Paper: instructions for use.

Cardboard: Kit box.

Plastic: test pack pouch, extraction buffer tube, cartridge

Several (not to be recycled): unused dipstick, unopened extraction buffer tube, swabs, dipstick pouch.

When used in accordance with Good Laboratory Practice, good occupational hygiene and instructions for use, the reagents should not present a hazard to health.

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1. Silver, Thomas Lee and Helayne. Etiology and Epidemiology of Preterm Premature Rupture of the Membranes. 2001. Clinics in Perinatology. Vol. 28, Issue 4.
2. Rutanen E-M, Pekonen F, Kärkkäinen T. Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clinica Chimica Acta (1993) 214: 73–81.
3. Rutanen E-M et al. Evaluation of a rapid strip test for insulin-like growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clin Chim Acta (1996) 253:91-101.
4. Ragosch V et al. Insulin like growth factor binding protein 1 (IGFBP-1) and fetales Fibronectin in der Diagnostik eines vorzeitigen Blasensprunges. GebFra (1996) 56:1-6.
5. Jain K and Morris P. A clinical study to evaluate the usefulness of the MAST test in diagnosing pre-labour rupture of membranes. J Obstet Gynaecol (1998) 18:33-36.
6. Kubota T and Takeuchi H. Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. J Obstet Gynecol Res (1998) 24:411-417.
7. Erdemoglu E and Mungan T. Significance of detecting insulin-like growth factor binding protein-1 in cervicovaginal secretions: Comparison with nitrazine test and amniotic fluid volume assessment. Acta Obstet Gynecol Scand (2004) 83:622-626.
8. H. Bohn and W. Kreaus. Isolierung und Charakterisierung eines neuen plazentaspezifischen Proteins (PP12). Arch. Gynecol (1980) 229: 279-291.

**Explanation of symbols • Vysvětlení symbolů • Symbolforklaring • Erklärung der Symbole • Explicación de los símbolos • Symbolien selitykset • Explication des symboles • Spiegazione dei simboli • Uitleg van symbolen • Symbolforklaringer • objaśnienie symboli • Explicação dos símbolos • Explicăția simbolurilor • Описание символов • Symbolförklaring**

**IVD**

In vitro diagnostic medical device  
Diagnostický zdravotnický prostředek in vitro  
Medicinsk udstyr til in vitro-diagnostik  
In-vitro-diagnostisches Medizinprodukt  
Producto sanitario para diagnóstico in vitro  
In vitro -diagnostiikkaan tarkoitettu lääkinällinen laite  
Appareil médical de diagnostic in vitro  
Dispositivo medico diagnóstico in vitro  
In-vitro diagnostisch medisch hulpmiddel  
In vitro-diagnostisk medisinsk enhet  
Wyrób medyczny do diagnostyki in vitro  
Dispositivo médico de diagnóstico in vitro  
Dispozitiv medical pentru diagnostic in vitro  
Медицинское изделие для диагностики in vitro  
Medicinteknisk produkt för in vitro-diagnostik

Consult instructions for use  
Přečtěte si návod k použití  
Se brugervejledningen  
Gebrauchsanweisung beachten  
Consulte las instrucciones de uso  
Katso käyttöohjetta  
Consulter le mode d'emploi  
Consultare le istruzioni per l'uso  
Raadpleeg de gebruiksaanwijzing  
Se bruksanvisningen  
Zajrzyj do instrukcji używania  
Consultar as instruções de utilização  
Consultați instrucțiunile de utilizare  
Обратитесь к инструкции по применению  
Se bruksanvisningen

Device for near-patient testing  
Prostředek pro vyšetření v blízkosti pacienta  
Enhed til patientnær test  
Produkt für patientennahe Tests  
Producto para diagnóstico analítico inmediato  
Vieritestauslaite  
Dispositif pour test au chevet du patient  
Dispositivo per analisi decentrate  
Medisch hulpmiddel voor NPT  
Enhet for pasientnær testing  
Wyrób do badań przytóżkowych  
Dispositivo para testes no ponto de tratamento  
Dispozitiv pentru testare în proximitatea pacientului  
Устройство для проведения анализов на месте лечения  
Produkt för patientnära testning

Sufficient for  
Postačuje pro  
Indeholder nok til  
Ausreichend für  
Suficiente evidencia clínica  
Lukumäärä  
Quantité suffisante pour  
Sufficiente per  
Toereikend voor  
Tilstrekkelig for  
Zawartość wystarczająca  
do <n> testów  
Suficiente para  
Suficient pentru  
Содержимого  
достаточно для  
проведения  
Tillräcklig för

Temperature limitation  
Teplotní omezení  
Temperaturgrænse  
Temperaturbegrenzung  
Limite de temperatura  
Lämpötilarajat  
Limite de température  
Limiti di temperatura  
Temperatuurbegrenzing  
Temperaturbegrensning  
Dopuszczalna  
temperatura  
Limitação de  
temperatura  
Limita de temperatură  
Температурный  
диапазон  
Temperaturbegrensning

Do not reuse  
Nepoužívejte opakovaně  
Må ikke genbruges  
Nicht wiederverwenden  
No reutilizar  
Kertakäyttöinen  
Ne pas réutiliser  
Non riutilizzare  
Niet opnieuw gebruiken  
Skal ikke gjenbrukes  
Nie używać powtórnie  
Não reutilizar  
A nu se reutiliza  
Запрет на повторное  
применение  
Får inte återvändas

**REF**

Catalogue number  
Katalogové číslo  
Katalognummer  
Katalognummer  
Número de catálogo  
Luettelonumero  
Número de catalogue  
Numero di catalogo  
Catalogusnummer  
Katalognummer  
Numer katalogowy  
Número de catálogo  
Număr de catalog  
Номер по каталогу  
Katalognummer

**LOT**

Batch code  
Kód šarže  
Partikode  
Chargencode  
Código de lote  
Eräkoodi  
Número de lot  
Codice lotto  
Partijcode  
Batchcode  
Kod partii  
Código do lote  
Număr de lot  
Код партии  
Parti-kod

Manufacturer  
Výrobce  
Producent  
Hersteller  
Fabricante  
Valmistaja  
Fabricant  
Produttore  
Fabrikant  
Producent  
Wytwórca  
Fabricante  
Producător  
Изготовитель  
Tillverkare

**IFU**

Instructions for use  
Návod k použití  
Brugervejledning  
Gebrauchsanweisung  
Instrucciones de uso  
Käyttöohje  
Mode d'emploi  
Istruzioni per l'uso  
Gebruiksaanwijzing  
Bruksanvisning  
Instrukcja użytkowania  
Instruções de utilização  
Instrucțiuni de utilizare  
Инструкции по  
применению  
Bruksanvisning

**SWAB**

Swab  
Tampón  
Podepind  
Tupfer  
Hisopo  
Näytteenottotikku  
Écouvillon  
Tampone  
Wattenstaafje  
Pinne  
Wymazówka  
Cotonete  
Tampón  
Тампон  
Provpinne

**EXTR**

Extraction buffer  
Extrakční pufr  
Ekstraktionsbuffer  
Extraktionspuffer  
Tampón de extracción  
Uuttopuskuri  
Tampon d'extraction  
Buffer di estrazione  
Extraktiebuffer  
Ekstraksjonsbuffer  
Bufor ekstrakcyjny  
Tampão de extração  
Soluție tampon de extracție  
Буфер для экстракции  
Extraktionsbuffert

**BUF**

**DIP**

Dipstick  
Měřicí tyčinka  
Extraktstrimmel  
Teststreifen  
Tira reactiva  
Testitikku  
Bandelette  
Striscia reagente  
Teststrip  
Teststrimmel  
Pasek testowy  
Vareta  
Bandă de testare  
Тест-полоска  
Teststicka

**CONT**

Contents  
Obsah  
Indhold  
Inhalt  
Contenido  
Sisältö  
Contenu  
Contenuto  
Inhoud  
Innhold  
Zawartość  
Conteúdo  
Conținut  
Содержимое  
Innehåll

**CART**

Cartridge  
Kazeta  
Patron  
Kartusche  
Cartucho  
Kasetti  
Cartouche  
Cartuccia  
Cartridge  
Kassett  
Wkład  
Cartucho  
Cartuş  
Картридж  
Kassett

Device not for self-testing  
 Prostredek není určen k sebetestování  
 Udstyret er ikke til selvtest  
 Produkt nicht zur Selbsttestung  
 Producto no destinado a autodiagnóstico  
 Ei itsetestaukseen  
 Dispositif non destiné à l'autotest  
 Dispositivo non per test autodiagnostici  
 Hulpmiddel niet voor zelftesten  
 Enhet ikke for selvtesting  
 Wyrób nieprzeznaczony do samokontroli  
 Dispositivo não destinado a autoteste  
 Dispozitiv care nu este destinat auto-testării  
 Устройство не для самостоятельного проведения анализов  
 Enheten är inte avsedd för självtestning

Keep away from rain  
 Chraňte před deštěm  
 Må ikke udsættes for regn  
 Vor Regen schützen  
 Mantener alejado de la lluvia  
 Suojaa sateelta  
 Tenir éloigné de la pluie  
 Tenere lontano dalla pioggia  
 Uit de buurt van regen houden  
 Skal ikke utsettes for regn  
 Chronić przed wilgocią  
 Manter afastado da chuva  
 A se feri de ploaie  
 Не допускать воздействия влаги  
 Håll borta från regn

Keep away from sunlight  
 Chraňte před slunečním zářením  
 Må ikke udsættes for sollys  
 Vor Sonnenlicht schützen  
 Mantener alejado de la luz solar  
 Suojaa suoralta auringonvalolta  
 Tenir éloigné des rayons du soleil  
 Tenere lontano dalla luce diretta del sole  
 Uit de buurt van zonlicht houden  
 Skal ikke utsettes for sollys  
 Trzymać z dala od światła słonecznego  
 Manter afastado da luz solar  
 A se feri de lumina soarelui  
 Не допускать попадания солнечного света  
 Håll borta från solljus

Single sterile barrier system  
 Systém jedné sterilní bariéry  
 Enkelt sterilt barriersystem  
 Einfaches Sterilbarriersystem  
 Sistema de barrera estéril único  
 Yksittäinen steriili suoja  
 Système à une barrière stérile  
 Sistema a barriera singola sterile  
 Enkelvoudig steriel barriersysteem  
 Enkelt, sterilt barriersystem  
 System pojedynczej bariery sterylnej  
 Sistema de barreira estéril única  
 Sistem de barieră sterilă unică  
 Система с одинарным стерильным барьером  
 Enkelt sterilt barriersystem

Authorised Representative in the European Community  
 Zplnomocněný zástupce v Evropském společenství  
 Autoriseret repræsentant i EU  
 Autorisierter Vertreter in der Europäischen Gemeinschaft  
 Representante autorizado en la Comunidad Europea  
 Euroopan yhteisössä toimiva valtettu edustaja  
 Représentant autorisé dans la Communauté européenne  
 Rappresentante autorizzato nella Comunità europea  
 Gemachtigd vertegenwoordiger in de Europese Gemeenschap  
 Autoriseret representant i EU  
 Autoryzowany przedstawiciel we Wspólnocie Europejskiej  
 Representante autorizado na Comunidade Europeia  
 Reprezentant autorizat în Comunitatea Europeană  
 Уполномоченный представитель в Европейском сообществе  
 Auktoriserad representant i Europeiska gemenskapen

Authorised representative in Switzerland  
 Zplnomocněný zástupce ve Švýcarsku  
 Autoriseret repræsentant i Schweiz  
 Autorisierter Vertreter in der Schweiz  
 Representante autorizado en Suiza  
 Valtuutettu edustaja Sveitsissä  
 Représentant autorisé en Suisse  
 Rappresentante autorizzato in Svizzera  
 Gemachtigd vertegenwoordiger in Zwitserland  
 Autoriseret representant i Sveits  
 Autoryzowany przedstawiciel w Szwajcarii  
 Representante autorizado na Suíça  
 Reprezentant autorizat în Elveția  
 Уполномоченный представитель в Швейцарии  
 Auktoriserad representant i Schweiz

Do not use if package is damaged  
 Nepoužívejte produkt, pokud má poškozený obal  
 Må ikke anvendes, hvis emballagen er beskadiget  
 Nicht verwenden, wenn Verpackung beschädigt ist  
 No utilizar si el envase está dañado  
 Käyttö kielletty jos pakkaus on vaurioitunut  
 Ne pas utiliser si l'emballage est endommagé  
 Non utilizzare se la confezione è danneggiata  
 Niet gebruiken als de verpakking beschadigd is  
 Skal ikke brukes hvis pakken er skadet  
 Nie używać, jeśli opakowanie jest uszkodzone  
 Não utilizar se a embalagem estiver danificada  
 A nu se utiliza dacă ambalajul este deteriorat  
 Не использовать при повреждении упаковки  
 Får inte användas om förpackningen är skadad

**MD**

Medical device  
 Zdravotnický prostředek  
 Medicinsk udstyr  
 Medizinprodukt  
 Producto sanitario  
 Lääkinnällinen laite  
 Appareil médical  
 Dispositivo medico  
 Medisch hulpmiddel  
 Medisinsk enhet  
 Wyrób medyczny  
 Dispositivo médico  
 Dispozitiv medical  
 Медицинское изделие  
 Medicinteknisk utrustning

PET

Polyethylene terephthalate  
 Polyethylentereftalát  
 Polyethylenterephthalat  
 Polyethylenterephthalat  
 Tereftalato de polietileno  
 Polyeteenitereftalaatti  
 Poly(téréphthalate d'éthylène)  
 Polietilene tereftalato  
 Polyethyleentereftalaat  
 Polyetylentereftalat  
 Poli(tereftalan etylenu)  
 Poli(tereftalato de etileno)  
 Polietilenä tereftalatä  
 Полиэтилентерефталат  
 Polyetentereftalat

**STERILEEO**

Sterilised using ethylene oxide  
 Sterilizováno pomocí ethylenoxidu  
 Steriliseret med ethylenoxid  
 Mit Ethylenoxid sterilisiert  
 Esterilizado con óxido de etileno  
 Steriloitu etyleenioksidilla  
 Stérilisé à l'oxyde d'éthylène  
 Sterilizzato con ossido di etilene  
 Gesteriliseerd met ethyleenoxide  
 Steriliseret med etylenoksid  
 Sterylizowany tlenkiem etylenu  
 Esterilizado com óxido de etileno  
 Sterilizat cu oxid de etilenă  
 Стерилизация оксидом этилена  
 Steriliserad med etylenoxid

PP

Polypropylene  
 Polypropylen  
 Polypropylen  
 Polypropylen  
 Polipropileno  
 Polypropeeni  
 Polypropylène  
 Polipropilene  
 Polypropyleen  
 Polypropylen  
 Polipropylen  
 Polipropileno  
 Polipropilena  
 Полипропилен  
 Polypropen

LDPE

Low density polyethylene  
 Polyethylen s nízkou hustotou  
 Polyethylen med lav densitet  
 Polyethylen niedriger Dichte  
 Polietileno de baja densidad  
 Matalatiheyksinen polyeteeni  
 Polyéthylène basse densité  
 Polietilene a bassa densità  
 Polyethyleen met lage dichtheid  
 Polyetylen med lav tetthet  
 Polietylen o niskiej gęstości  
 Polietileno de baixa densidade  
 Polietilenä de joasä densitate  
 Полиэтилен низкой плотности  
 Polyeten med låg densitet

PAP

Paper  
 Papír  
 Papir  
 Papier  
 Papir  
 Papel  
 Paperi  
 Papier  
 Carta  
 Papier  
 Papier  
 Papier  
 Hârtie  
 Бумага  
 Papper

PAP

Non-corrugated fibreboard  
 Nevlnitá dřevovláknitá deska  
 Pap uden bølger  
 Nicht gewellte Pappe  
 Cartón no ondulado  
 Kartonki  
 Panneau de fibres non ondulé  
 Pannello di fibra non ondulato  
 Vezelplaat zonder golfkarton  
 Glatt papp  
 Płyta pilśniowa niekarbowana  
 Painel de fibras não ondulado  
 Panouri neondulate  
 Негофрированный картон  
 Kartong (ej wellpapp)

Warning  
 Varování  
 Advarsel  
 Warnhinweis  
 Advertencia  
 Varoitus  
 Avertissement  
 Attenzione  
 Waarschuwing  
 Advarsel  
 Ostrzeżenie  
 Advertência  
 Avertisment  
 Осторожно  
 Varning

CONT MCI/MI

Contains <0,01 % reaction mass of: 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) [EC No: 611-341-5]  
Obsahuje < 0,01 % hmotnosti reakčních látek: 5-chlor-2-methyl-2H-isothiazol-3-on a 2-methyl-2H-isothiazol-3-on (v poměru 3 : 1) [č. EC: 611-341-5]  
Inneholder <0,01 % reaktionsmasse af: 5-chlor-2-methyl-2H-isothiazol-3-on og 2-methyl-2H-isothiazol-3-on (3:1) [EF-nr.: 611-341-5]  
Enthält <0,01 % Reaktionsmasse von: 5-Chlor-2-methyl-2H-isothiazolin-3-on und 2-Methyl-2H-isothiazolin-3-on (3:1) [EG-Nr.: 611-341-5]  
Contiene <0,01 % de masa reactiva de: 5-cloro-2-metil-2H-isotiazol-3-ona y 2-metil-2H-isotiazol-3-ona (3:1) [n.º de EC: 611-341-5]  
Sisältää 5-Kloori-2-metyyli-2H-isotiatol-3-ionin [EY nro. 247-500-7] ja 2-metyyli-2H-isotiatol-3-ionin [EY nro. 220-239-6] (3:1) seosta  
Contient une masse de réaction < 0,01 % de : 5-chloro-2-méthyl-2H-isothiazol-3-one et de 2-méthyl-2H-isothiazol-3-one (3:1) [N° CE : 611-341-5]  
Contiene <0,01 % della massa di reazione di: 5-cloro-2-metil-2H-isotiazol-3-ona e 2-metil-2H-isotiazol-3-ona (3:1) [CE N.: 611-341-5]  
Bevat < 0,01 % reactiemassa van: 5-chloor-2-methyl-2H-isothiazol-3-on en 2-methyl-2H-isothiazol-3-on (3:1) [EC-nr: 611-341-5]  
Inneholder < 0,01 % reaksjonsmasse av: 5-kloro-2-metyl-2H-isotiazol-3-one og 2-metyl-2H-isotiazol-3-one (3:1) [EC-nr.: 611-341-5]  
Zawiera <0,01% masy reakcyjnej: 5-chloro-2-metylo-2H-izotiazol-3-onu i 2-metylo-2H-izotiazol-3-onu (w stosunku 3:1) [nr WE: 611-341-5]  
Contém massa de reação <0,01% de: 5-cloro-2-metil-2H-isotiazol-3-ona e 2-metil-2H-isotiazol-3-ona (3:1) [N.º CE: 611-341-5]  
Conține <0,01% masă de reacție alcătuită din: 5-cloro-2-metil-2H-izotiazol-3-onă și 2-metil-2H-izotiazol-3-onă (3:1) [nr. CE: 611-341-5]  
Содержит < 0,01 % реакционной массы 5-хлор-2-метил-2H-изотиазол-3-она и 2-метил-2H-изотиазол-3-она (3 : 1) [номер EC: 611-341-5]  
Innehåller <0,01 % reaktionsmassa av: 5-kloro-2-metyl-2H-isothiazol-3-on och 2-metyl-2H-isothiazol-3-on (3:1) [EG-nr: 611-341-5]

This product fulfills the requirements of the European Parliament and of the Council on *in vitro* diagnostic medical devices  
Tento produkt splňuje požadavky Evropského parlamentu a Rady na diagnostické zdravotnické prostředky *in vitro*  
Dette produkt opfylder kravene fra Europa-Parlamentet og Rådet om medicinsk udstyr til *in vitro*-diagnostik  
Dieses Produkt erfüllt die Anforderungen des Europäischen Parlaments und des Rates zu Medizinprodukten für die In-Vitro-Diagnose  
Este producto cumple con los requisitos del Parlamento Europeo y del Consejo sobre productos sanitarios de diagnóstico *in vitro*  
Tämä tuote täyttää Euroopan parlamentin ja neuvoston *in vitro*-diagnostiikkaan tarkoitettujen lääkinnällisten laitteiden vaatimukset  
Ce produit est conforme aux exigences du Parlement européen et du Conseil des appareils médicaux de diagnostic *in vitro*  
Questo prodotto soddisfa i requisiti del Parlamento europeo e del Consiglio europeo in materia di dispositivi medico-diagnostici *in vitro*  
Dit product voldoet aan de vereisten van het Europees Parlement en de Raad betreffende *in vitro* diagnostische medische hulpmiddelen  
Dette produktet oppfyller kravene fra Europaparlamentet og Rådet til *in vitro*-diagnostiske medisinske enheter  
Ten produkt spełnia wymagania Parlamentu Europejskiego i Rady dotyczące wyrobów medycznych do diagnostyki *in vitro*  
Este produto cumpre os requisitos do Parlamento Europeu e do Conselho para dispositivos médicos de diagnóstico *in vitro*  
Acest produs îndeplinește condițiile stabilite de Parlamentul European și de Consiliul UE privind dispozitivele medicale pentru diagnosticul *in vitro*  
Этот продукт соответствует требованиям Европейского парламента и Совета в отношении медицинских изделий для диагностики *in vitro*  
Den här produkten uppfyller kraven i Europaparlamentets och rådets direktiv om medicintekniska produkter för *in vitro*-diagnostik

SWAB

Puritan Medical Products Company  
31 School Street, Guilford, Maine 04443-0149, USA

EC REP

EMERGO EUROPE Arnhem, The Netherlands

Actim<sup>®</sup>PROM  
1ngeni

Actim Oy  
Klovinpellontie 3, FI-02180 Espoo, Finland  
Tel. +358 9 547 680  
[www.actimtest.com](http://www.actimtest.com)  
[actim@actimtest.com](mailto:actim@actimtest.com)

© Actim Oy. All rights reserved.  
Actim<sup>®</sup> is a trademark of Actim Oy.

AORACE30831-2 2024-05-02