

# Actim<sup>®</sup> Calprotectin

A convenient rapid test for  
the diagnosis and monitoring  
of intestinal inflammation.



# Actim<sup>®</sup> Calprotectin

## Technical specifications

Test principle	Semiquantitative immuno-chromatographic dipstick test
Target protein	Calprotectin
Specimen type	Feces
Reading time	10 minutes
Interpretation	Visual
Detection limits	< 50 µg/g 50-200 µg/g >200 µg/g
Test components	1 dipstick in a perforating device. Specimen Dilution Buffer, 3 ml
Number of tests/kit	20 tests (35032 ETAC)
Storage	+2...+30 °C
Sample storage	Original feces sample can be stored 4 days in refrigerator at +2...+8 °C. The specimen suspension can be stored for one week at +2...+25 °C.



## Actim Calprotectin Controls kit (35000 ETAC) for quality control

Kit components	1 negative control (35080 ETAC) 1 positive control, low (35081 ETAC) 1 positive control, high (35082 ETAC) Reconstitution Solution (30889 ETAC)
Storage	+2...+8 °C

## Ordering information

Product	Product code
Actim Calprotectin, 20 tests	35032 ETAC
Actim Calprotectin Controls kit	35000 ETAC

**Actim Oy**  
– part of Medix Biochemica Group  
Klovinpellontie 3, FI-02180 Espoo, Finland  
actim@actimtest.com  
[www.actimtest.com](http://www.actimtest.com)

