

Actim[®] PROM

A highly specific test for
the detection of premature
rupture of fetal membranes.



Actim[®] PROM

Technical specifications

Test principle	Qualitative immuno-chromatographic test
Target protein	IGFBP-1 (insulin-like growth factor binding protein-1)
Specimen type	Vaginal secretions
Specimen collection	Swab specimen from posterior fornix with or without a sterile speculum
Reading time	5 minutes
Interpretation	Visual
Cut-off	25 µg/l in the extracted specimen
Test components	1 swab, 1 tube of Specimen Extraction Solution, 1 dipstick
Number of tests/kit	20 tests (30832ETAC) 10 tests (30831ETAC) 1 test (30830ETAC)
Storage	+2... +25°C



Actim PROM Controls kit (30800ETAC) for quality control

Kit components	1 negative control (30880ETAC) 1 low positive control (30881ETAC) 1 high positive control (30882ETAC) 1 reconstitution solution (30889ETAC)
Storage	+2... +25°C

Ordering information

Product	Product code
Actim PROM, 20 tests	30832ETAC
Actim PROM, 10 tests	30831ETAC
Actim PROM, 1 test	30830ETAC
Actim PROM Controls	30800ETAC

Actim Oy

– part of Medix Biochemica Group
Klovinpellontie 3, FI-02180 Espoo, Finland
actim@actimtest.com

www.actimtest.com

