

Actim[®] Partus

A rapid diagnostic test
to help assess the risk
of preterm birth.



Actim[®] Partus

Technical specifications

Test principle	Qualitative immuno-chromatographic test
Target protein	phIGFBP-1 (phosphorylated form of insulin-like growth factor binding protein-1)
Specimen type	Cervical secretions
Specimen collection	Swab specimen from cervical os. during speculum examination
Reading time	5 minutes
Interpretation	Visual
Cut-off	10 µg/l in the extracted specimen
Test components	1 swab, 1 tube of extraction buffer , 1 dipstick
Number of tests/ kit	10 tests (31931ETAC) 1 test (31930ETAC)
Storage	+2... +25°C



Actim Partus Controls kit (31900ETAC) for quality control

Kit components	1 negative control (31980ETAC) 1 low positive control (31981ETAC) 1 high positive control (31982ETAC) 1 reconstitution solution (30889ETAC)
Storage	+2... +25°C

Product information

Product	Product code
Actim Partus, 10 tests	31931ETAC
Actim Partus, 1 test	31930ETAC
Actim Partus Controls	31900ETAC

Actim Oy

– part of Medix Biochemica Group
Klovinpellontie 3, FI-02180 Espoo, Finland
actim@actimtest.com

www.actimtest.com

